

## The Ultimate Inside Lift Solution

Fully-powered to hoist and transport a power chair or scooter inside a cargo area. Triple axis lift provides fully-powered lifting, boom extension and rotation for the easiest, most effortless loading.

### Features & Benefits

- 12 full inches of power extension to reach beyond rear bumper
- Powered lifting, extension and rotation
- New, simplified base fits in a smaller footprint
- Easier to install and disassemble
- Faster mounting with fewer holes, less drilling when necessary
- Passenger or driver side install
- Wireless key fobs included
- 3-year transferable warranty



Specifications	AL435 Axis III Lift
Accommodates	Power chairs, scooters, manual wheelchairs
Lifting Capacity	400 lbs. (181 kg)
Installed Weight	76 lbs. (34.5 kg)
Adjustable Boom Height	30" – 41"
Adjustable Boom Length	32.5" – 43.5"
Adjustable Boom Pitch	2° - 35°
Securement	Docking device included
Warranty	3-year transferable

### Key Features



Convenient wireless key fobs available



Easy-to-use hand control with Soft-Touch buttons.