Fits Virtually All Power Chairs

Harmar's AL500 and AL500 Heavy Duty accommodates virtually any power chair. Available in traditional 350 lbs. or heavy-duty 400 lbs. models.

Features & Benefits


• Easy, one-switch operation

Self-tensioning Q'Straint retractors; no chair modifications required

• 27.5" x 39" aluminum deck

Manual crank backup and license plate mount included

Swing-Away option available



800.833.0478

| Specifications | AL500 | AL500 HD Heavy Duty |
|------------------|--|--|
| Accommodates | Power chairs | Power chairs |
| Platform | 27.5" x 39"aluminum deck | 27.5" x 39"aluminum deck |
| Lifting Capacity | 350 lbs. (160 kg) | 400 lbs (181 kg) |
| Installed Weight | 75 lbs. (34 kg) | 77 lbs. (35 kg) |
| Securement | Self-tensioning Q'Straint retractors | Self-tensioning Q'Straint retractors |
| Hitch | Class II or greater | Heavy Duty Class III only |
| Options | Battery pack, Swing-Away, Power chair covers | Battery pack, Swing-Away, Power chair covers |
| Warranty | 3-year transferable | 3-year transferable |

