



2014 Truck Lifts



				
AL835 Pickup 350 HD	AL825 Pickup 225 Speed	AL825CC Crew Cab 225 HD	AL815CC Crew Cab 150 Speed	AL435T Tailgater

Accommodates	Standard and Extended Cab Trucks	Standard and Extended Cab Trucks	Crew Cab Trucks	Crew Cab Trucks	All Pick-up Trucks
Mobility Devices	Manual Chairs, Scooters and Small Power Chairs up to 350 lbs.	Manual Chairs, Scooters and Small Power Chairs up to 225 lbs.	Manual Chairs, Scooters and Small Power Chairs up to 225 lbs.	Manual Chairs, Scooters and Small Power Chairs up to 135 lbs.	All Scooters and Power Chairs up to 400 Lbs.
Adjustable Boom Height	34" - 41"	34" - 41"	34" - 41"	34" - 41"	31" - 39"
Adjustable Boom Length	28.75" - 46.5"	28.75" - 46.5"	46.5" - 67"	46.5" - 67"	21" - 33"
Adjustable Boom Pitch	N/A	N/A	N/A	N/A	0° - 35°
Lifting Weight Capacity	350 lbs. (159 kg)	225 lbs. (102 kg)	225 lbs. (102 kg)	135 lbs. (61 kg)	400 lbs. (181 kg)
Installed Weights	116 lbs. (53 kg)	116 lbs. (53 kg)	121 lbs. (55 kg)	121 lbs. (55 kg)	81 lbs. (37 kg)
Docking Device	Included	Included	Included	Included	Included

Benefits



MADE IN USA



2075 47th Street,
Sarasota, FL 34234
800-833-0478
www.harmar.com

- | | | | | |
|---|---|--|--|---|
| <ul style="list-style-type: none"> • Heavy Duty Model for Standard / Extended Cab Trucks • Powered Lifting & Rotation • Offset Base Mounting Fits Tightly to Sidewall and Fender walls • Driver or Passenger Side Mounting Bases Available • Dashboard Mounted On/Off Switch • Remote Hand Control (Wireless Option) • Dual-strapped Hook Design Prevents Mobility Device from Spinning During Load/Unload • 3-Year Transferable Warranty | <ul style="list-style-type: none"> • Twice the Speed of AL835 for Standard / Extended Cab Trucks • Powered Lifting & Rotation • Offset Base Mounting Fits Tightly to Sidewall and Fender walls • Driver or Passenger Side Mounting Bases Available • Dashboard Mounted On/Off Switch • Remote Hand Control (Wireless Option) • Dual-strapped Hook Design Prevents Mobility Device from Spinning During Load/Unload • 3-Year Transferable Warranty | <ul style="list-style-type: none"> • Crew Cab Model • Powered Lifting & Rotation • Offset Base Mounting Fits Tightly to Sidewall and Fender walls • Driver or Passenger Side Mounting Bases Available • Dashboard Mounted On/Off Switch • Remote Hand Control (Wireless Option) • Dual-strapped Hook Design Prevents Mobility Device from Spinning During Load/Unload • 3-Year Transferable Warranty | <ul style="list-style-type: none"> • Twice the Speed of AL825CC for Crew Cab Trucks • Powered Lifting & Rotation • Offset Base Mounting Fits Tightly to Sidewall and Fender walls • Driver or Passenger Side Mounting Bases Available • Dashboard Mounted On/Off Switch • Remote Hand Control (Wireless Option) • Dual-strapped Hook Design Prevents Mobility Device from Spinning During Load/Unload • 3-Year Transferable Warranty | <ul style="list-style-type: none"> • Unique Offset Arm Loads/Unloads Mobility Device from Curbside or over 24" Tailgate • Power Extension Reaches 43" Over Tailgate • Passenger-side Install Only • Wireless Remote Control • 3-Year Transferable Warranty |
|---|---|--|--|---|