



## Pickup 225 Speed Lift

Fits regular or extended cab trucks for manual wheelchairs, small power chairs, and scooters. Speed model is two times faster than heavy-duty models and ideal for manual wheelchairs.

Specifications	AL825
Accommodates	Power chairs & scooters
Lifting Capacity	210 lbs. (95 kg)
Installed Weight	116 lbs. (53 kg)
Powered	Lifting & rotation
Adjustable Boom Height	34" – 41"
Adjustable Boom Length	28.75" – 46.5"
Securement	Docking device included
Available Options	Battery pack, wireless remote
Warranty	3-year transferable

## Features & Benefits

- *Delivers chair close to driver or passenger side door*
- *Powered lifting and rotation*
- *Offset base mounting fits tightly to sidewall & fender walls*
- *Only two motors to control with remote hand control*
- *Driver or passenger-side mounting bases available*
- *Dashboard-mounted on/off switch*
- *Dual-strapped hook prevents mobility device from spinning during load/unload*
- *Remote hand control included; wireless option available*
- *Docking device included*

