

AL585SEL & AL585SS

SPECIALLY-DESIGNED LIFTS FOR FRONT-WHEEL-DRIVE POWER CHAIRS

AL585SEL FWD Power Chair Lift

Designed
for FWD
power chairs
(Driver's side
entry only)



AL585SS FWD Power Chair Lift

Designed
for FWD
power chairs
(Passenger's side
entry only)



- Drive-on, drive-off platform
- Platform automatically folds when not in use
- Lightweight: Just 76 lbs. installed
- Keyed, one switch operation
- Requires no chair modifications
- Manual crank backup
- Padded hold-down arm safely and automatically holds chair

SPECIFICATIONS

Lifting Capacity:	350 lbs. (160 kg) Upgradable to 400 lbs.
Aluminum Deck Size:	27.5" x 39"
Installed Lift Weight:	75 lbs.
Hitch Height (range):	5" – 18"
Hitch Class:	II or III
Warranty:	3-year transferable

