



Harmar® SL350OD Outdoor Stair Lift

Withstands the Elements

Harmar's Outdoor Stair Lift is engineered with a durable extruded aluminum track and marine-grade vinyl seating to handle rain, wind, snow, and temperatures from 0 to 125 degrees. Includes header, internal chassis, and large weatherproof covers for even more protection.



Features & Benefits

- Designed for long-term outdoor use
- Constant-pressure rocker controls
- Operates on standard 115 VAC household current
- Seat and footrest fold to save space
- Seat cushions of exterior marine-grade vinyl
- 350 lbs. capacity



Specifications	SL350OD
Maximum Weight Capacity:	350 lbs. (160 kg)
Overspeed Governor:	No
Track (rail) Type:	Extruded Aluminum
Travel:	16' Standard; 25' Maximum (7.6 m)
Average Number of Return Trips per Charge: (varies with load, length)	N/A
Symmetrical Lift: (Ships and installs non handed)	No
Control in Armrest: (left or right hand operation)	Yes
Need to Open Chassis to Install or Change Sides:	Yes
Lift Mounts to Steps or Wall:	Steps
Minimum Folded Width:	14" (35.5 cm)
Minimum Footrest Height:	4" (10 cm)
Clear Distance Between Armrests:	19" (48.3 cm)
Floor to Seat Height Low / High:	20" - 24" (50.8 cm - 61 cm)
Minimum Wall to Stair Side of Rail:	10" (25.4 cm)
Seat Depth:	14" (35.6 cm)
Backrest Height From Top of Seat:	14" (35.5 cm)
Electrical Requirements:	120 VAC 15A
Operation Power:	AC
Speed:	20 fpm
Incline Range:	30 - 45°
Drive System:	Cable Drive
Warranty:	2-Years on Components
Safety Devices:	Limit Switches Final Limit Switch Slack Cable Brake w/ Switch Footrest Obstruction Switch Swivel Seat Cut-off Switch Constant Pressure Controls Seat Belt

Available Seat Colors:



Light Almond (Standard)



Pine Green



Mocha



Harmar® 800-833-0478 | www.harmar.com
America's Lift Leader™

