

## Characteristics

Tested for use with 12 common chemotherapy drugs, with all drug permeation data printed on the box. This medium weight nitrile glove is extremely elastic and comfortable while providing the tensile strength required for today's demanding procedures. Polymer inner surface for effortless donning. Textured, slightly tacky finish for excellent gripping ability even in the presence of moisture. Aqua blue color.



**Exam Glove  
Non-Sterile**

**NitriDerm®  
AquaBLUE**

**Nitrile  
Series 197**



## PRODUCT DETAILS

SIZE	ITEM NO.	PACKAGING	DESCRIPTION
XS	197050	100 Gloves/box, 10 boxes/case	Gloves, Exam, Nitrile, Non Sterile, Powder-Free, Textured
S	197100	100 Gloves/box, 10 boxes/case	
M	197200	100 Gloves/box, 10 boxes/case	
L	197300	100 Gloves/box, 10 boxes/case	
XL	197350	100 Gloves/box, 10 boxes/case	
XXL	197400	90 Gloves/box, 10 boxes/case	

## Features

- Low Modulus
- Non-Latex
- Textured Finish
- CHEMO Rated

## Benefits

- Softer, More Comfortable Fit
- No Risk of Latex Allergens
- Improved Wet/Dry Grip
- Tested for CHEMO Use

*Product Solutions You Trust*