

# INOGEN ONE<sup>®</sup>G2 System



Inogen One<sup>®</sup> G2  
AC Power Supply\*  
Item # BA-301



Inogen One<sup>®</sup> G2  
DC Power Supply\*  
Item # BA-302



Inogen One<sup>®</sup> G2  
Cart & Carry Bag\*  
Item # CA200



Inogen One<sup>®</sup> G2  
12-Cell Lithium Ion Battery\*  
Item #BA200



Inogen One<sup>®</sup> G2  
24-Cell Lithium Ion Battery\*\*  
Item #BA224



## The Inogen One<sup>®</sup> G2 System (IS-200):

- Only 7.0 Pounds
- 6 Flow Setting
- Extremely Quiet
- Authorized by FAA onboard aircraft

\*Included with system

\*\*Accessory sold separately

## Oxygen. Anytime. Anywhere®

The new Inogen One® G2 delivers the independence of our original Inogen One portable concentrator, but it's smaller, lighter, and makes more oxygen with a longer battery life. It's a portable concentrator small enough to hold in your hand. That means there are no tanks to refill and you never have to worry about running out of oxygen.

The Inogen One® G2 is your one solution, your single solution, for oxygen at home, away, and for travel (approved for use on all US airlines). The Inogen One® G2 is Oxygen. Anytime. Anywhere®.

**To learn more, call toll-free 1-866-543-4664.**

<b>Weight</b>	7.0 lbs. (includes 12-cell battery)
<b>Size</b>	Length: 10.7" Width: 3.9" Height: 9.5" (including handle)
<b>Oxygen Flow</b>	Pulse dose delivery system Six flow settings from 1-6
<b>Power Sensing</b>	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power: for mobile use in car/airplane
<b>Battery</b>	Duration (12-cell battery): Up to 5 hours Recharge approximately 4 hours with AC or DC power Duration (24-cell battery): Up to 10 hours Recharge approximately 8 hours with AC or DC power
<b>Noise Level</b>	38 Decibels*
<b>Warranty</b>	3 year or lifetime limited warranty**
<b>Operation</b>	Simple control functions and easy-to-read LCD display
<b>Use</b>	Designed for stationary, portable, and travel

\*At setting 2

\*\*See Inogen One® G2 Limited Warranty Statement for complete warranty information

Manufactured by:

**inogen®**

326 Bollay Drive  
Goleta, California 93117  
1-866-543-4664  
[www.inogenone.com](http://www.inogenone.com)  
[info@inogen.net](mailto:info@inogen.net)

