

INOGEN ONE^{G3} System



Inogen One[®]
AC Power Supply*
Item # BA-301



Inogen One[®] G3
DC Power Cable*
Item # BA-306



Inogen One[®] G3
Carry Bag*
Item # CA-300



**Only
4.8 pounds!**



Inogen One[®] G3
Single Lithium Ion Battery*
Up to 4.5 hours run time
Item #BA-300



Inogen One[®] G3
Double Lithium Ion Battery**
Up to 9 hours run time
Item #BA-316

The Inogen One[®] G3 System (IS-300):

- 45% Smaller
- 30% Lighter
- 27% More Oxygen Per Pound
Compared to the Inogen One G2

*Included with system

**Accessory sold separately

Oxygen. Anytime. Anywhere[®]

The new Inogen One[®] G3 delivers the independence of previous Inogen One portable concentrators, but it's smaller, lighter, and makes more oxygen per pound. It's a portable concentrator light enough to be worn care-free by the active oxygen user. That means there are no tanks to refill and you never have to worry about running out of oxygen.

The Inogen One[®] G3 is your one solution, your single solution, for oxygen at home or away. The Inogen One[®] G3 is Oxygen. Anytime. Anywhere[®].

To learn more, call toll-free 1-800-461-6751.

| | |
|----------------------|--|
| Weight | 4.8 lbs. (includes single battery) |
| Size | Length: 8.75" Width: 3.0" Height: 7.25" |
| Oxygen Flow | Pulse dose delivery system Four flow settings 1, 2, 3, 4 |
| Power Sensing | AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car |
| Battery | Duration (single battery): Up to 4.5 hours Recharge approximately 2-4 hours with AC or DC power Duration (double battery): Up to 9 hours Recharge approximately 4-8 hours with AC or DC power |
| Noise Level | 42 Decibels* |
| Warranty | 3 year or lifetime limited warranty** |
| Operation | Simple control functions and easy-to-read LCD display |
| Use | Designed for 24/7 use |

*At setting 2

**See Inogen One[®] Limited Warranty Statement for complete warranty information

Manufactured by:

 inogen[®]

326 Bollay Drive
Goleta, California 93117
1-800-461-6751
www.inogenone.com
info@inogen.net

