

NEW
FROM **Inova**Labs 

LifeChoice **Activox™ 4L**



Lightweight at **4.8 lbs.**

Up to **4 LPMeq**
PULSE-WAVE™ Delivery

Advanced **O₂ Sensor**
standard on every unit


Up to **10.25 hours** of
internal battery runtime

FAA-Approved

Compatible with
Activox DUO2™



**Tank
Free!**



Enhanced Transflective Display

- > Clear alarm & charging notifications
- > Simplified hour meter & purity readings
- > Tailored compliance tracking

> STANDARD PACKAGE



Accessories:

4-Way Carry Case
Adjustable Straps
AC Power Supply
DC Power Supply
Accessory Bag
Nasal Cannula

External
Battery
Available*



> KEY SPECIFICATIONS

Weight	4.8 lbs. (2 kg)			
Oxygen Delivery	PULSE-WAVE™ Delivery			
Oxygen Concentration	90% (+/-3%)			
Oxygen Concentration Sensor Alarm	85% or less			
Pulse Settings	1 LPMeq	2 LPMeq	3 LPMeq	4 LPMeq
Internal Battery Runtime*	10.25 Hrs	8.25 Hrs	5 Hrs	4 Hrs
External Battery Runtime*†	4.75 Hrs	3.75 Hrs	2.75 Hrs	2.25 Hrs
Charge Time*	4 Hrs Internal, 2 Hrs External			
Average Sound Level at 2 LPMeq	44 dB at 3 ft. (1 m)			

*Hours are approximate. Runtimes are based on breath per minute.

†Not included in Standard Package.

Manufactured & Distributed By InovaLabs 

1.800.220.0977 | www.InovaLabs.com

3500 Comsouth Drive, Suite 100, Austin, TX 78744 USA

© 2014 Inova Labs, Inc. All rights reserved. Made in the USA. I2M-405 Rev A Oct 2014