


INVACARE®
VERANDA™ 4000
LIGHTWEIGHT WHEELCHAIR



Where clinical efficacy meets practicality. This elegantly designed lightweight wheelchair provides comfort and mobility for the patient in an economical package. Complete with center-mount front riggings packaged with the chair and a space-saving design, the **Invacare® Veranda™ 4000 Wheelchair** is a practical solution for the clinical provider.

Features

Center mount pivot footrests. Footrests swing inward under the chair, providing space around the chair for transfer. Simplicity of inventory and ease-of-use for the patient.

Center mount elevating articulating legrests. The elevating articulating legrests pivot and extend at the knee providing a more comfortable and natural articulating motion.

Tool-free adjustability. Tool-free height adjustable back with three back heights, and tool-free footrest length adjustment, provide a simple and versatile solution.

Cross-compatibility. Interchangeability of service parts helps cut inventory and increase cash flow.

Smaller package design. The chair comes with the back off and is assembled without tools allowing for a smaller footprint in your warehouse and delivery van.

HCPCS Code: K0004



1. Center mount elevating articulating legrest
2. Tool-less height adjustable back
3. Natural motion legrests
4. Center mount pivot footrests
5. Footrests swing inward under the chair
6. Seating options available



Available in 16", 18" and 20" seat widths, and 16" and 18" seat depths with footrests or legrests in the box.

Model #	
4V66FFR	16x16 Frame w/ Desk Length Flipback Arms and Footrest
4V66FLR	16x16 Frame w/ Desk Length Flipback Arms and Elevating Articulating Legrests
4V86FFR	18x16 Frame w/ Desk Length Flipback Arms and Footrest
4V86FLR	18x16 Frame w/ Desk Length Flipback Arms and Elevating Articulating Legrests
4V88FFR	18x18 Frame w/ Desk Length Flipback Arms and Footrest
4V88FLR	18x18 Frame w/ Desk Length Flipback Arms and Elevating Articulating Legrests
4V06FFR	20x16 Frame w/ Desk Length Flipback Arms and Footrest
4V06FLR	20x16 Frame w/ Desk Length Flipback Arms and Elevating Articulating Legrests
4V08FFR	20x18 Frame w/ Desk Length Flipback Arms and Footrest

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2013

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form No. 13-553
130874



Model #	Overall Width	Overall Height	Overall Depth w/o riggings	Folded Width	Seat Width	Seat Depth	Seat-to-Floor Height
V16, V18, V20	23.5", 25.5", 27.5"	33.5"-37.6"	31.5" - 33.2"	14"	16", 18", 20"	16", 18"	17.5", 19.5"
Back Height	Caster Height	Wheel Height	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty	
V16, V18, V20	16", 17", 18"	8"	24"	300 lb.	54 lb.	33"x12"x 36"	1 Year