Invacare® Bariatric Drop-Arm Commode Features

Model # 6599

1000 lb. weight capacity

Easy-to-operate drop-arm mechanism

Extra-wide snap-on seat

Removable padded backrest

Adjustable height legs

Durable steel frame design

Comes complete with pail and lid

HCPCS Code: E0168



1. Bariatric size



Invacare® Heavy-Duty Commode Features

Model # 6497

450 lb. weight capacity

Designed for bedside use

Rectangular pan and pan holder

Armrests for added comfort

HCPCS Code: E0168



1. Bariatric size









- 1. Invacare® Pan and Holder
 Model #13245X000 (pan)
 - Model #13245X000 (pan) Model # 1027004 (holder) 8-quart pan capacity
- 2. Invacare® Pail with Lid

Model # 6317 (12 per carton) 12-quart pail capacity

3. Invacare® Splash Shield

Model #6319 (12 per carton)















Model #	Dimensions HxWxD E	Width Between Arms	Seat Surface Dimensions	Seat Opening Dimensions	Seat-to-Floor Height	Seat-to-Arm Height	Weight
9630-1, 9630-4	23"-29"x25"x19"	19"	16"Wx14"D	8.5"Wx10"D	16"-22"	7.5"	12 lb.
9650-4	23"-29"x25"x19"	19"	16"Wx14"D	8.5"Wx10"D	16"-22"	7.5"	7 lb.
9610-4	23"-29"x25"x19"	19"	16"Wx14"D	8.5"Wx10"D	16"-22"	7.5"	12 lb.
9669	27"-31.5"x20.25"x17.	5" 17.5"	14"Wx15"D	8.25"Wx9.5"D	17"-21"	8"	22 lb.
9830-4	25"-31"x21.5"x18.5"	18"	13.25"Wx14.5"[8.5"Wx9.75"D	16.5"-21.5"	8"	13 lb.
6500-BHD	25"-29"x36"x21"	30"	27.75"Wx18.5"[9"Wx12"D	17"-21"	8.25"	20 lb.
6599	27"-31"x28"x20.5"	25"	23.5"Wx18"D	8.5"Wx12"D	19.25"-23.25"	7.75"	21 lb.
6497	25"-31"x30"x19"	21 75"	15 5"Wx15 25"Γ) 8 25W"x12"D	16 25"-22"	9"	16 lb







Weight	Units	Limited
Capacity	Per Carton	Warranty
300 lb.	4	5 Years
300 lb.	4	5 Years
1,000 lb.	1	5 Years
300 lb.	4	5 Years
250 lb.	1	5 Years
650 lb.	2	5 Years
300 lb.	4	5 Years
450 lb.	2	5 Years
	Capacity 300 lb. 300 lb. 1,000 lb. 300 lb. 250 lb. 650 lb. 300 lb.	Capacity Per Carton 300 lb. 4 300 lb. 4 1,000 lb. 1 300 lb. 4 250 lb. 1 650 lb. 2 300 lb. 4

Invacare Corporation

www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

©2012

Invacare Corporation.

All rights reserved.

Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Specifications are subject to change without notification.

Form No. 07-056

Rev. 03/12

111276



INVACARE® COMMODES





Invacare Commodes offer consumers the comfort and stability they need. The frame is lightweight and consumers gain additional comfort and support from the 16" wide seat. The Invacare All-in-One Commode can be used bedside or can act as a toilet safety frame or raised toilet seat.

Invacare® I•Class™ All-In-One Commode Features

Model # 9630-1, 9630-4

Can be used bedside or with backrest removed, can act as a toilet safety frame or raised toilet seat

Large 16" wide snap-on seat with lid for added comfort and support

Comes complete with commode bucket, bucket cover and splash shield

Wing nuts allow for removal of backrest when needed

HCPCS Code: E0163



1. 16" wide snap-on seat for added comfort and support



Invacare® All-In-One Aluminum Commode Features

Model # 9650-4

Lightweight aluminum frame is rust resistant

Can be used bedside or, with backrest removed, can act as a toilet safety frame or raised toilet seat

Large 16" wide snap-on seat with lid for added comfort and support

Comes complete with commode bucket, bucket cover and splash shield

Wing nuts allow for removal of backrest when needed

HCPCS Code: E0163



1. 16" wide snap-on seat for added comfort and support



Invacare® Safeguard Commode Features

Model # 9610-4

Gray-coated steel

Can be used as a raised toilet seat or a toilet safety frame

Equipped with a splash shield

Snap-on seat

HCPCS Code: E0163



Invacare® Drop-Arm Commode Features

Model # 9669

Can be used over the toilet with optional splash shield

Accommodates transfers in restricted areas

Drop-arm lever operates simply: lift to lock or push to release

Gray powder coated, 1" welded steel tubing

Snap-on seat and lid

HCPCS Code: E0163



Invacare® Folding Commode Features

Model # 9830-4

Offering all the great features that you expect. It folds up in seconds for storage or to transport and comes with armrests, angled legs for stability, non-marring rubber tips, anti-rattle collars that reduce noise, snap-on removable back, pail, splash shield and lid.

HCPCS Code: E0165



Invacare® Bariatric Commode Features

Model # 6500-BHD

650 lb. weight capacity

Extra-wide, snap-on seat

Rectangular pan can be removed from the front or rear of the commode

Adjustable seat-to-floor height

HCPCS Code: E0168



1. Bariatric size

