

INVACARE® ASPIRATOR



The **Invacare® Aspirator** is a durable suction unit designed for most home care suction requirements.

Features

High-performance vacuum specifications

Adjustable vacuums from 0" to 23" Hg

Durable piston vacuum pump

Easy-to-read gauge

Non-marring suction cup feet

Single-use collection jars and connecting tubing for quick and easy replacement (1 set included)

UL and CSA listed

HCPCS Code: E0600



Suction cup feet

Disposable collection jar, Model # IRC1140 (HCPCS Code: A7002)

Tube and filter kit (not shown) Model # 1148804

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2012

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
®. All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject
to change without
notification.
Form No. 11-252
111290



All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.