

INVACARE®
PRONTO® M41™
POWER WHEELCHAIR



The **Invacare® Pronto® M41™ Power Wheelchair** features the latest technology in electronics and in-line motors. This combination provides increased efficiency and battery range, as well as an increased top speed. The M41 Wheelchair also features a center-wheel drive platform, simplified Invacare® SureStep® Suspension and six wheels on the ground providing stability and maneuverability similar to other Pronto models. All of these features are packed into a simple, clean looking and easy-to-service compact frame that provides mobility with style for people on the go.

Features

Compact size for a tight turning radius
 High efficiency in-line motors for greater speed and longer battery range
 Charger connection conveniently located on joystick
 Stylish design
 3-amp off-board charger

Three seat options:

Medium-back fold-down seat (18" W x 17" D)

High-back semi-recline office-style seat

High-back semi-recline captain's seat with solid seat pan

Auto style positioning strap

HCPCS Code: K0822, K0823



1. Center-Wheel Drive
 2. Simplified Invacare® SureStep® Suspension



Invacare Corporation
 www.invacare.com

USA

One Invacare Way
 Elyria, Ohio
 44035-4190
 (800) 333-6900

Canada

570 Matheson Blvd. E.,
 Unit 8
 Mississauga, Ontario
 L4Z 4G4 Canada
 (800) 668-5324

©2012

Invacare Corporation.
 All rights reserved.
 Trademarks are identified
 by the symbols ™ and
 ®. All trademarks are
 owned by or licensed
 to Invacare Corporation
 unless otherwise noted.
 Specifications are subject
 to change without
 notification.
 Form No. 06-033
 Rev. 09/12
 120347



Model #	Overall Dimensions	Seat Dimensions	Seat-to-Floor Height	Drive Wheels	Ground Clearance	Turning Radius	Incline
M41FDR/B, M41RSOLID16R/B, M41RSOLIDR/B, M41RSOLID20R/B, M41SR16R/B M41SRR/B M41SR20R/B	33.7" L 24.2" W 33.1 H	16"-20" W 16"-20" D	19"-23" (18"-22" w/o cushion)	10.5"	2"	17.75" Front 19.5" Rear	7°
	Total Weight	Weight Capacity	Motor	Maximum Speed	Battery	Battery Range	Limited Warranty
	157 lb.	300 lb.	2 Pole	4.8 mph	U1	Up To 14.5 Miles	5 Years Base Frame 3 Years Seat Frame 1 Year Electronics/ Motor/Gearbox

Do not operate on roads, streets or highways.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.