# INVACARE® PRONTO® M94™ POWER WHEELCHAIR





The Invacare® Pronto® M94™ Power Wheelchair with Invacare® SureStep® Suspension combines innovative features and exceptional quality to provide a power wheelchair that meets the mobility needs of bariatric consumers. The M94 Wheelchair features the patented Invacare® SureStep® Suspension that makes negotiating thresholds and transitions of up to 3" high effortless. Also featured is the patented Stability Lock technology that helps ensure that the chair will remain stable with the consumer in control. These two innovative features paired with its sleek styling make the Pronto M94 Wheelchair a truly exceptional power chair for consumers weighing up to 500 lb.

## **Captain's Seat Features**

Semi-reclining back (standard) for additional comfort

Width adjustable flip-back arms (standard) with optional height adjustment

Ergonomic seating contours provide support for longer durations of use

Adjustable height headrest for head and neck support Built-in seat pouch on back

Standard easy to clean grey vinyl upholstery 18", 20", 22" and 24" widths

#### **Rehab Seat Features**

Angle adjustable back offers range of 80° to 100° Sling back upholstery provides comfort and support for mild positioning needs

Round canes accept other back options for additional positioning

## HCPCS Code: K0852, K0853











- 1. True Center-Wheel Drive
- 2. Invacare® SureStep® Suspension
- 3. Stability Lock
- 4. Invacare® TCD™ Traction Control Design
- 5. Invacare® MK6i™ Electronics



# **Invacare Corporation**

www.invacare.com

### USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

# Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

©2012

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols <sup>TM</sup> and (B). All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
Specifications are subject to change without notification.
Form No. 04-001
Rev. 09/12
120351















	190P	904	900	( <del>0</del> 0₽	906	-=-	( de la composición della comp
Model #	Overall Dimensions	Seat Dimensions	Seat-to-Floor Height	Drive Wheels	Ground Clearance	Turning Radius	Incline
M94	43" L* 31.5" W 51" H**	18"-24" W 16"-22" D	19"-24"	14"	3.25"	21.5"	9°



Total

Weight

290 lb



Weight

Capacity

500 lb.



Motor

4-Pole HD



Maximum

Speed

5.1 mph



Battery

22NF



Battery

Range

Up To

11.95 Miles



Limited

Warranty

5 Years



\*\* For other seating option specifications refer to the Invacare Visual Configurator Do not operate on roads, streets or highways.

Base Frame/ 3 Years Seat Frame 1 Year Electronics/ Motor/Gearbox



All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.