



INVACARE®
PRONTO® 31
POWER WHEELCHAIR



The **Invacare® Pronto® 31 Power Wheelchair** continues the rich legacy of the Pronto series. This compact front-wheel drive power wheelchair features simple handling, durability and comfort all in an ultra-sleek design. Two new refreshing colors, Blue Mist and Black Fog make the Pronto 31 Wheelchair the model of modern day power mobility.



Features

Comfortable 18" captain's seat with 360° rotation and 45° of adjustability
Seat and back covers are easily replaceable
Durable, molded-in-color shroud available in Blue Mist and Black Fog
Compact size for tight turning, ideal for in-home use
Charger connection conveniently located on joystick

Simple, intuitive user controls
Refurbishment Kits allow for economical refreshing
13 month extended warranty on motors, electronics and select components

HCPCS Code: K0816

















- 1. Tool-less removable shroud
- 2. Removable, full length, flip-up, adjustable arms
- 3. Adjustable angle and height fold-down footboard
- 4. 9" flat free drive wheels
- 5. 50 Amp PG VR2 controller



Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Model #	 Overall Dimensions	 Seat Dimensions	 Seat-to-Floor Height	 Drive Wheels	 Ground Clearance	 Turning Radius	 Incline
P31BLUE, P31BLACK	38.75" L* 23.5" W	18" W 18" D	20.5" - 22"	9"	1.75"	26.5"	6°
	 Total Weight	 Weight Capacity	 Motor	 Maximum Speed	 Battery	 Battery Range	 Limited Warranty
	153 lb.	300 lb.	2 Pole	4 mph	U1	Up To 13 Miles	5 Years Base Frame 3 Years Seat Frame 13 Months Motors and Electronics 6 Months Batteries

*Footboard down
Do not operate on roads, streets or highways.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

©2014
Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®.
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form No. 13-572
Rev. 03/14
140260

