


INVACARE®  
**VERANDA™**  
WHEELCHAIR

---



---

The **Invacare® Veranda™ Wheelchair** provides comfort and mobility in an economical package. Complete with front riggings packaged with the chair, the Invacare Veranda wheelchair combines style and durability at an affordable price. Available in 16", 18" and 20" with footrests and/or legrests.

---

## Features

Stylish powder-coated frame includes front riggings

- Footrests match chair
- Legrests are chrome

Carry pocket on back


Large step-tube end cap

Padded nylon upholstery

Comfortable padded arms

Dual axle position

**HCPCS Code: K0001/K0002**



**Ask For It!**



1. Front riggings
2. Nylon upholstery
3. Dual axle position
4. Seating options available

**Invacare Corporation**  
www.invacare.com

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

©2012

Invacare Corporation.  
All rights reserved.  
Trademarks are identified  
by the symbols ™ and  
®. All trademarks are  
owned by or licensed  
to Invacare Corporation  
unless otherwise noted.  
Specifications are subject  
to change without  
notification.  
Form No. 07-064  
Rev. 03/12  
111238



Model #	Overall Width	Overall Height	Overall Depth w/o riggings	Folded Width	Seat Width	Seat Depth	Seat-to-Floor Height
V16	23"	36"	30"	14"	16"	16"	17.5", 19.5"
V18	25"	36"	30"	14"	18"	16"	17.5", 19.5"
V20	27"	36"	30"	14"	20"	16"	17.5", 19.5"
	Back Height	Caster Height	Wheel Height	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty
V16	18"	8"	24"	250 lb.	54 lb.	33"x12"x 36"	1 Year
V18	18"	8"	24"	250 lb.	54 lb.	33"x12"x 36"	1 Year
V20	18"	8"	24"	250 lb.	54 lb.	33"x12"x 36"	1 Year

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.