



INVACARE® TOP END®  
**EXCELERATOR™ SERIES**  
HANDCYCLE



If you want a great way to exercise, cross-train or just have fun, one of the **Invacare® Top End® Excelerator™ Series Handcycles** is just what you need. This upright and extremely maneuverable three-wheel handcycle is capable of speeds up to 15 mph. The smooth seven-speed Shimano® hub with its reversing drum brake delivers hands-on braking and allows you to tackle hills or level terrain with equal ease. The oval tube carbon steel frame construction with adjustable sliding seat makes transfers in and out easier too.

The Excelerator Handcycle is built to your specifications, plus you can order the popular mountain drive option which adds 7 extra low gears for climbing hills or the V crankset option with its wider cranks to allow more clearance for those who are tall or require additional space. No excuses...Get Cranking!

---

## Features

Carbon steel frame and fork

Adjustable fore/aft seat

Adjustable back height

Fold over back upholstery

Adjustable crank height

Chain tension idler

Full chain guard, reflectors and safety flag

Adjustable height/angle footrest system

Shimano® 7 speed hub w/reversing drum brake

Parking brake

175 mm crankset w/one chain ring

Ovalized ergonomic hand pedals

Performance wheels w/cruiser tires

Quick release axles


9° camber

Cushion and hook and loop seat restraint

Weight capacity 350 lb.

---

## Model #: XCL



1. Twist Shifter for 7 speed hub

2. Quad Twist shifter option

3. Parking brake

4. V crankset option

5. Mountain drive transmission option

6. Frictionless ovalized handpedal option



---

## Options

Mountain drive transmission adds 7 low gears

V crank allows maximum clearance

Quad twist shifter helps those with hand impairments

Tri-pin quad handpedals

Crank width adapter set

Quad gloves

Water bottle and cage

Tow bar for wheelchair

Cordless computer

Front safety light

Rear flashing light

Helmet

Seat restraint: Auto style fastener

Mirror

Crutch holder and strap

Bike rack for trailer hitch

Top End t-shirt

---

## Features

Carbon steel frame and fork  
14" adjustable seat and back height  
Fold over back upholstery  
Adjustable crank height  
Chain tension idler  
Full chain guard, reflectors and safety flag  
Adjustable height/angle footrest system  
Shimano® 7 speed hub with reversing drum brake

Parking brake  
155 mm crankset with one chain ring  
Foam covered horizontal handpedals  
20" Cruiser wheels and tires  
Quick release axles  
9° camber  
Cushion and hook and loop seat restraint  
Weight capacity 250 lb.

---

## Model #: LXC



1. Parking brake
2. Horizontal handpedals
3. Hook and loop safety seat restraint

---

## Options


Quad twist shifter option  
Quad gloves  
Mirror  
Water bottle and cage  
Helmet  
Front safety light  
Rear flashing light  
Crutch holder and strap  
Bike rack for trailer hitch  
Seat restraint: Auto style fastener

---




### 1 V Crankset

The Invacare® Top End® V crankset option accommodates consumers who might have "clearance" issues. Not compatible with mountain drive option



### 2 Mountain Drive Transmission

Available on the Excelerator handcycle  
 Adds 7 low gears (7 low + 7 standard = 14 total)  
 Easy shifting, especially for quadriplegics – just push-in the center button  
 Lightweight planetary transmission  
 A must have for hilly terrain



### 3 Full Chain Guard

A full chain guard on all 7 speed handcycles is a great safety feature

**Invacare Corporation**  
[www.invacare.com](http://www.invacare.com)

#### USA

One Invacare Way  
 Elyria, Ohio  
 44035-4190  
 (800) 333-6900

#### Canada

570 Matheson Blvd. E.,  
 Unit 8  
 Mississauga, Ontario  
 L4Z 4G4 Canada  
 (800) 668-5324

#### Top End

[www.topendwheelchair.com](http://www.topendwheelchair.com)  
 Sports and Recreation  
 Products  
 4501 63rd Circle North  
 Pinellas Park, Florida  
 33781-5914  
 (727) 522-8677  
 (800) 532-8677  
 Fax: (727) 522-1007

©2012

Invacare Corporation.  
 All rights reserved.  
 Trademarks are identified  
 by the symbols ™ and ®.  
 All trademarks are owned  
 by or licensed to Invacare  
 Corporation unless  
 otherwise noted. Shimano  
 is a registered trademark of  
 Shimano Inc.  
 Form No. 12-509  
 120736

