

VIP-515 Transport Class

TRANSPORTER

TILT IN SPACE CHAIR

SPECIFICATIONS	VIP-515-TP16	VIP-515-TP18
Seat Width	16 in.	18 in.
Seat Depth	17 in.	17 in.
Seat Height	20 in.	20 in.
Back Height	18+ in.	18+ in.
Overall Height	53 in.	53 in.
Overall Open Width	24 in.	26 in.
Weight w/o Riggings	33 lbs.	33 lbs.
Weight Capacity	250 lbs.	250 lbs.
HCPCS Code	E1038	E1038

(SUBJECT TO CHANGE WITHOUT NOTICE)



FIND THE **PERFECT** FIT FOR YOU...

- Only lightweight tilt-in-space chair that also folds!
- Shifts pressure to reduce bedsores
- Intelligent 0-35 center-of-gravity weight shifting to shift pressure and reduce bedsores
- Center-Of-Gravity weight shifting
- Specially designed ergonomic head pillow is included
- Paint coating uses special Anti-bacterial materials
- With newly designed swing in and away footrests
- Flip-back armrests that are height adjustable flip-back height-adjustable armrests
- Anti-bacterial and anti-staining Aegis upholstery Weight
- Full-length adjustable height removable armrests
- Seat Belt Restraint Included