

## **ADULT CO<sub>2</sub> SAMPLING CANNULA**

Soft and flexible to ensure your patient's comfort


# ADULT CO<sub>2</sub> SAMPLING CANNULA

## For Simultaneous CO<sub>2</sub> Monitoring & Oxygen Delivery

- CO<sub>2</sub> sampling is collected through nasal prongs for accurate monitoring.
- Oxygen is delivered through the bottom ports, below the nasal prongs.
- Cannula is fitted with a universal connector designed to fit any flowmeter without the need for adapters.
- Both male, female and universal luer connectors are available.



Simultaneous CO<sub>2</sub> Monitoring & Oxygen Delivery



Universal O<sub>2</sub> connector



Male luer connector



Female luer connector



Universal luer connector



Item #	Size	Oxygen Tubing	O <sub>2</sub> Connector	Gas Sampling Line	Luer Connector	Adaptor	Pkg.
HCS4560	Adult	7 ft.	Universal	4" pigtail	Male	22mm ID	25 ea/cs
HCS4561	Adult	7 ft.	Universal	4" pigtail	Female	22mm ID	25 ea/cs
HCS4568	Adult	7 ft.	Universal	10 ft.	Male	22mm ID	25 ea/cs
HCS4569	Adult	7 ft.	Universal	10 ft.	Female	22mm ID	25 ea/cs
HCS4565	Adult	10 ft.	Universal	10 ft.	Universal	22mm ID	25 ea/cs
HCS4566	Adult	14 ft.	Universal	14 ft.	Universal	22mm ID	25 ea/cs
VTB4107	Adult	7 ft.	Standard	7 ft.	Universal	22mm ID	25 ea/cs
VTB4107F	Adult	7 ft.	Universal	7 ft.	Universal	22mm ID	25 ea/cs
HCS4562	Pediatric	7 ft.	Universal	4" pigtail	Male	22mm ID	25 ea/cs
HCS4563	Pediatric	7 ft.	Universal	4" pigtail	Female	22mm ID	25 ea/cs



**Medline Industries, Inc.**  
One Medline Place  
Mundelein, IL 60060

**Medline United States**  
**1-800-MEDLINE (633-5463)**

www.medline.com | info@medline.com

**Medline Canada**  
**1-800-396-6996**

www.medline.ca | canada@medline.com

**Medline México**  
**01-800-831-0898**

www.medlinemexico.com | mexico@medline.com

Some products may not be available for sale in Mexico or Canada. Please contact your account representative for additional information.  
We reserve the right to correct any errors that may occur on this sell sheet. ©2011 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc.