

GUARDIAN® Performance and Style Rollators

- A superior group of rollators for a customer looking for more than the basic
- Engineered with design features that solve specific problems for the user
- Unique colors and additional accessories make these models ideal for a retail showroom
- All models use loop-style brakes you push down to lock and pull up to release

New Colors!



MDS86826SC
SERIES

New Colors!



MDS86825SLF*



MDS808200TR

MDS86826 SERIES

Deluxe Comfort Rollator with Memory Foam Seat

Find comfortable seating and smooth operation in one model. Extra thick, memory foam seat contours comfortably to create an ideal cushion for tired users while the large, 8" (20 cm) wheels glide smoothly over most surfaces. Convenient cup holder and cane holder come standard. Also comes standard with front storage pouch!

MDS86825SL SERIES

Freedom Ultra Light Rollator

Weighing only 11 lbs. (5 kg), the Freedom rollator is the lightest model on the market. Light weight and folding capabilities make this the ideal model to get in and out of the car/storage. Height-adjustable arms and legs let users find a perfect fit as this rollator accommodates users 4'11"-6'4" (150-193 cm).

MDS808200TR SERIES

Excel Translator: Combination Transport Chair and Rollator

In just seconds, our Excel Translator easily converts from a rollator to a transport chair and vice-versa. Comes with a convenient cup holder, under-the-seat storage bag, zippered side carrying case, and easy-to-adjust push-button foot rests. You will be delighted with the comfortable, spacious seat which is more than 75% roomier than an average rollator seat.

Item Number	Name	Color	Wt. Cap. Lbs (kg)	Wheels Inch (cm)	Inch (cm) Between Handles	Seat Height Min./Max. Inch (cm)	Height Min./Max. Inch (cm)	Product Weight Lbs (kg)	Pkg
MDS86826SC	Deluxe Comfort	Midnight Blue	300 (136)	8 20)	17 (43)	23 (58)	32.5-37.5 (83-95)	17 (8)	Each r
MDS86826SCG*	Deluxe Comfort	Green	300 (136)	8 20)	17 (43)	23 (58)	32.5-37.5 (83-95)	17 (8)	Each r
MDS86826SCR*	Deluxe Comfort	Red	300 (136)	8 20)	17 (43)	23 (58)	32.5-37.5 (83-95)	17 (8)	Each r
MDS86825SLR	Freedom	Burgundy	250 (113)	6 15)	17 (43)	18-23 (46-58)	29-36 (74-91)	11 (5)	Each r
MDS86825SLB	Freedom	Smoky Blue	250 (113)	6 15)	17 (43)	18-23 (46-58)	29-36 (74-91)	11 (5)	Each r
MDS86825SLBLK	Freedom	Black	250 (113)	6 15)	17 (43)	18-23 (46-58)	29-36 (74-91)	11 (5)	Each r
MDS86825SLT	Freedom	Titanium	250 (113)	6 15)	17 (43)	18-23 (46-58)	29-36 (74-91)	11 (5)	Each r
MDS86825SLWT	Freedom	Pearl White	250 (113)	6 15)	17 (43)	18-23 (46-58)	29-36 (74-91)	11 (5)	Each r
MDS86825SLF*	Freedom	Fushia	250 (113)	6 15)	17 (43)	18-23 (46-58)	29-36 (74-91)	11 (5)	Each r
MDS86825SLG*	Freedom	Lime Green	250 (113)	6 15)	17 (43)	18-23 (46-58)	29-36 (74-91)	11 (5)	Each r
MDS86825SLP*	Freedom	Salmon	250 (113)	6 15)	17 (43)	18-23 (46-58)	29-36 (74-91)	11 (5)	Each r
MDS808200TRR	Translator	Red	250 (113)	8 20)	18 (46)	19.5 (50)	32-36.5 (81-93)	18 (8)	Each
MDS808200TR	Translator	Blue	250 (113)	8 20)	18 (46)	19 (48)	32-36.5 (81-93)	18 (8)	Each

*Coming Soon!