



244 Bariatric Treatment Table



Barrier-Free® Accessibility - Safe and Efficient Care for your Patient, you and your Staff

850-pound Weight Capacity - With an 850-pound patient weight capacity, the Ritter 244 Bariatric Power Treatment Table is ready, willing and able to safely lift special needs patients every time.

18-inch Low Height - With an incredible Barrier-Free® 18-inch low height, the Ritter 244 is designed to provide a solution to the all too common question, "How are we going to safely and effectively get our special needs patients onto and off of the treatment table?"

By lowering the seat to 18 inches from the floor, the Ritter 244 performs this sometimes difficult task. It becomes a very simple routine, requiring minimal assistance from you and your staff. The patient's dignity is maintained and the likelihood of injury is greatly reduced, providing efficient, effective patient care.

34-inch High Height - The Ritter 244 raises the seat to a height of 34 inches, designed to bring the patient to you while helping to minimize the stress and strain on your back. By having the patient in the right position, the exam can be completed properly and efficiently.

The Ritter 244 features a 32-inch wide seat that offers unsurpassed patient comfort and support for your special needs patients. The seamless design of the upholstery looks great and offers an easy-to-clean surface. The thick foam padding with built-in pillow provides the comfort your patients deserve.



With an ergonomic and intuitive hand control, patient positioning is simple and efficient.





The optional adjustable arm system accessory is designed to offer a sense of security for patients and a stable, secure platform to assist them in getting onto and off of the table. With the simple touch of a lever, the arms can be rotated and locked into one of three positions so the arms are always in the right place.





18" Low Height 34" High Height



*Product shown with optional adjustable arm system accessory.

PECIFICATIONS





Fully Powered, **Multi-Position Movement**

In addition to fully powered height adjustment and powered back and tilt, angle adjustments are designed to allow for a thorough examination and safe positioning. Safety switches in the base section stop the table from moving if a slight amount of force is applied to the base cover.

Ritter 244 Bariatric Power Treatment Table

An optional foot control accessory is designed to provide a sterile field and an alternate method to the hand control for positioning patients.



Height: Back Support Angle:

$Maximum = 15^{\circ}$ Upholstery Dimensions: Seat and Back Area: 32" W (81.3 cm)

 $Minimum = -15^{\circ}$

Specifications

Patient Weight Capacity:

Minimum 18" (45.7 cm)

Maximum 34" (86.4 cm

850 lb (385 kg)

 $Minimum = 0^{\circ}$

Seat Tilt Angle:

 $Maximum = 65^{\circ}$

Head and Foot Area: 30" W (76.2 cm)

Overall Length: 80" (203.2 cm)

Floor Clearance from Floor to Bottom of Base:

Leveling Feet Installed: 2.75" (7 cm)

Caster Accessory Installed: 3.5" (8.9 cm)

Paper Roll Holder Capacity: 18", 21" or 24" (45.7 cm, 53.3 cm or 61.0 cm

Accessories

Adjustable arm system (right and left side arms included)

Casters

Foot control

3' power cord

Electrical

120 VAC, 2.5 Amps, 60 Hz

Safety Standards UL 60601-1 and CAN/CSA C22.2 No. 601.1 - M90

Ritter Upholstery Colors:

Navy-231
Dusty Blue-233
Moss-230
Clay-234
Pebble Grey-216
Shadow-232
Perfect Plum-229
Black-312

CARB 93120.2 Phase 2 Compliant

Midmark is an ISO 13485 and ISO 9001 Certified Company.

For more information or a demonstration, contact your Midmark Dealer or Midmark Corporation, Versailles, Ohio. Call: 1-800-MIDMARK Fax: 1-800-365-8631 midmark.com



The optional casters permit easy transportation of your table. The casters can support a 650 lb stationary patient load but should not be used to move patients.



© Midmark Corporation 2004 Manufactured for Midmark Corporation. Products subject to improvement changes without notice. Litho in U.S.A. 007-0537-00 Rev. G1 (5/16)