



Reston self-adhering foam products are designed to protect against skin damage from splints, casts, and prostheses. It may be used as padding for such devices as beds, chairs, and operating room tables.

- Versatile padding adheres well to a variety of surfaces.
- Lightweight and conformable padding can be cut, layered and shaped, allowing it to be adapted for many applications.
- Available in medium- and high-support pads or light support foam roll to meet your application needs.

3M™ Reston™ Self-Adhering Foam Products



Product Description:

3M™ Reston™ Self-Adhering Foam Products consist of a porous polyurethane foam covered on one side with a hypoallergenic, acrylate adhesive and a protective, easy to remove, silicone coated paper release liner.

Contraindications:

This product is not to be used as a dressing in post-surgical treatment procedures such as liposuction.

Warnings:

Reston self-adhering foam products have not been evaluated for use as a compression dressing material in surgical procedures such as liposuction. The product should not be applied to the skin under tension. Such uses can result in damage to the skin such as blistering or tearing which may lead to hypopigmentation and the risk of serious infection. Reston self-adhering foam products should not be placed in direct contact with an open wound.

Ordering Information:

Catalog No.	Product Information	Size	Pads/Package	Packages/Case
3M™ Reston™ Self-Adhering Foam Pad				
 1560M	Medium Support	7/16 in thick 7 7/8 in x 11 3/4 in 20 cm x 29,8 cm	10	5
 1561H	High Support	15/16 in thick 7 7/8 in x 11 3/4 in 20 cm x 29,8 cm	5	5
Catalog No.	Product Information	Size	Rolls/Package	Packages/Case
3M™ Reston™ Self-Adhering Foam Roll				
 1563L	Light Support	3/16 in thick 4 in x 196 in 10 cm x 5 cm	1	5

If you have any questions or comments, in the USA please contact your local Skin & Wound Care Specialist, the 3M Health Care Customer Helpline at **1-800-228-3957** or visit **www.3M.com/skinhealth**. Outside the United States, contact your local 3M subsidiary.



Skin & Wound Care Division
3M Health Care
 3M Center, Building 275-4W-02
 St. Paul, MN 55144-1000
 USA
 1-800-228-3957
www.3M.com/healthcare

3M Canada Company
 P.O. Box 5757
 London, Ontario N6A 4T1
 Canada
 1-800-364-3577
www.3M.com/ca/healthcare

3M and Reston
 are trademarks of 3M.
 Please recycle.
 Printed in U.S.A.
 Used under license in Canada.
 © 3M 2009. All rights reserved.
 70-2008-0838-7