


3M™ Tegaderm™ +Pad Film Dressing with Non-Adherent Pad 3M™ Medipore™ +Pad Soft Cloth Adhesive Wound Dressing



- All-in-one sterile dressing convenience
- Highly conformable and comfortable
- Unique non-adherent pad



Trusted Performance Delivered Two Styles

3M™ Tegaderm™ +Pad Film Dressing with Non-Adherent Pad

A waterproof, absorbent dressing.

- All-in-one sterile dressing—reduces application time
- A bacterial and viral* barrier—helps prevent external bacteria and contaminants from entering the site
- A waterproof barrier—patient may bathe or shower without removing dressing
- Conforms to body contours, providing comfort
- Unique, non-adherent absorbent pad

**In vitro* testing shows that the transparent film of Tegaderm™, Tegaderm HP™ and Tegaderm™ CHG dressings provides a viral barrier from viruses 27 nm in diameter or larger while the dressing remains intact without leakage.



Tegaderm™ +Pad on epidural



Tegaderm™ +Pad on hip incision



Tegaderm™ +Pad on arthroscopic incision



Tegaderm™ +Pad
*Full range of sizes to fit
every application need.*

Suggested Applications

Primary Dressing Applications

- Post-surgical – especially on areas where swelling or movement is expected
- Acute wounds such as large abrasions, lacerations, and burns
- Light-to-moderately-draining wounds
- Superficial and partial thickness wounds
- I.V. catheter sites
- Fragile skin

in

of Absorbent Island Dressings

3M™ Medipore™ +Pad Soft Cloth Adhesive Wound Dressing

A gentle, conformable, absorbent dressing.

- All-in-one sterile dressing—reduces application time
- Conforms to body contours, providing comfort
- Gentle, porous and breathable
- Unique, non-adherent absorbent pad



Medipore™ +Pad on IV catheter site



Medipore™ +Pad on knee



Medipore™ +Pad on appendectomy

Secondary Dressing Applications

- Over hydrogels, alginates, gauze, or other primary dressings



Medipore™ +Pad
Full range of sizes to fit
every application need.

Catalog No. NDC No.	Dressing Size	Pad Size	Dressings/ Box	Boxes/ Case	HCPCS* Code
3M™ Tegaderm™ +Pad Film Dressing with Non-Adherent Pad					
3582 NDC 8333-3582-01	2 in. x 2-3/4 in. 5 cm x 7 cm	1 in. x 1-1/2 in. 2,5 cm x 4 cm	50	4	A6203
3584 NDC 8333-3584-01	2-3/8 in. x 4 in. 6 cm x 10 cm	1 in. x 2-3/8 in. 2,5 cm x 6 cm	50	4	A6203
3586 NDC 8333-3586-01	3-1/2 in. x 4 in. 9 cm x 10 cm	1-3/4 in. x 2-3/8 in. 4,5 cm x 6 cm	25	4	A6203
3587 Oval NDC 8333-3587-01	3-1/2 in. x 4-1/8 in. 9 cm x 10,5 cm	1-3/4 in. x 2-3/8 in. 4,5 cm x 6 cm	25	4	A6203
3588 NDC 8333-3588-01	6 in. x 6 in. 15 cm x 15 cm	4 in. x 4 in. 10 cm x 10 cm	25	4	A6203
3589 NDC 8333-3589-01	3-1/2 in. x 6 in. 9 cm x 15 cm	1-3/4 in. x 4 in. 4,5 cm x 10 cm	25	4	A6203
3590 NDC 8333-3590-01	3-1/2 in. x 8 in. 9 cm x 20 cm	1-3/4 in. x 6 in. 4,5 cm x 15 cm	25	4	A6203
3591 NDC 8333-3591-01	3-1/2 in. x 10 in. 9 cm x 25 cm	1-3/4 in. x 8 in. 4,5 cm x 20 cm	25	4	A6203
3593 NDC 8333-3593-01	3-1/2 in. x 13-3/4 in. 9 cm x 35 cm	1-3/4 in. x 11-3/4 in. 4,5 cm x 30 cm	25	4	A6204

Catalog No. NDC No.	Dressing Size	Pad Size	Dressings/ Box	Boxes/ Case	HCPCS* Code
3M™ Medipore™ +Pad Soft Cloth Adhesive Wound Dressing					
3562 NDC 8333-3562-01	2 in. x 2-3/4 in. 5 cm x 7 cm	1 in. x 1-1/2 in. 2,5 cm x 4 cm	50	4	A6219
3564 NDC 8333-3564-01	2-3/8 in. x 4 in. 6 cm x 10 cm	1 in. x 2-3/8 in. 2,5 cm x 6 cm	50	4	A6219
3566 NDC 8333-3566-01	3-1/2 in. x 4 in. 6 cm x 10 cm	1-3/4 in. x 2-3/8 in. 4,5 cm x 6 cm	25	4	A6219
3568 NDC 8333-3568-01	6 in. x 6 in. 15 cm x 15 cm	4 in. x 4 in. 10 cm x 10 cm	25	4	A6219
3569 NDC 8333-3569-01	3-1/2 in. x 6 in. 9 cm x 15 cm	1-3/4 in. x 4 in. 4,5 cm x 10 cm	25	4	A6219
3570 NDC 8333-3570-01	3-1/2 in. x 8 in. 9 cm x 20 cm	1-3/4 in. x 6 in. 4,5 cm x 15 cm	25	4	A6219
3571 NDC 8333-3571-01	3-1/2 in. x 10 in. 9 cm x 25 cm	1-3/4 in. x 8 in. 4,5 cm x 20 cm	25	4	A6219
3573 NDC 8333-3573-01	3-1/2 in. x 13-3/4 in. 9 cm x 35 cm	1-3/4 in. x 11-3/4 in. 4,5 cm x 30 cm	25	4	A6220

* HCPCS codes as referenced in the Local Coverage Determination Surgical Dressing Policy.

Disclaimer: HCPCS codes have been provided to assist you in the preparation of Medicare Part B claims. Please note, however, that the reimbursement information provided by 3M Health Care and its representatives is intended to provide general information relevant to coverage and coding for 3M products. Insurers' reimbursement policies can vary and the use of the codes discussed here does not guarantee that an insurer will cover or pay at any particular level. Health care providers should exercise independent clinical judgment in choosing the codes which most accurately describe the products provided.



Skin & Wound Care Division
3M Health Care
 3M Center, Building 275-4W-02
 St. Paul, MN 55144-1000
 USA
 1-800-228-3957
www.3M.com/healthcare

3M Canada
 Post Office Box 5757
 London, Ontario N6A 4T1
 Canada
 1-800-364-3577

3M, Medipore and Tegaderm
 are trademarks of 3M.
 Please recycle.
 Printed in U.S.A.
 © 3M 2008, 2010. All rights reserved.
 70-2008-8925-4 (Bolger)