

Steel Transport Chair

The Everest & Jennings Transport Chair is designed as an affordable patient transfer device. Its compact size and lighter weight allow a companion to lift and store the chair with ease.

Features:

- Available in Silver Vein Hammertone finish
- Fold down back makes transport and storage easy
- Black, padded, Nylon upholstery is flame retardant
- 17" and 19" seat widths
- Carry straps on seat allow for ease of lifting and folding
- Footrests are industry compatible and offer easy tool free adjustment
- Full length padded arms provide user comfort
- Weighs 24 lb w/o front rigging
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of Medicare code E1038

Warranty:

1 year on frame and crossbraces



| Item Number | Description | Seat | Seat Depth | Seat Height | Back Height | Overall Depth | Overall Height | Overall Width | Folded Width | Weight |
|-------------|---------------------|------|------------|-------------|-------------|---------------|----------------|---------------|--------------|--------|
| EJ795-1 | Steel - Silver Vein | 17" | 15 3/4" | 19" | 18" | 37" | 39" | 24" | 10" | 24 lb |
| EJ796-1 | Steel - Silver Vein | 19" | 15 3/4" | 19" | 18" | 37" | 39" | 24" | 10" | 24 lb |

NOTE: Everest & Jennings maintains the policy of continual product improvement, therefore, we reserve the right to make changes without notice. All dimensions are ±1/4". *Without Front Rigging.

