

Sprague

Stethoscopes




Sprague Stethoscope

The Pro Advantage sprague stethoscope features a double-tube configuration for improved stereo separation as well as a patented chestpiece valve stem which delivers outstanding acoustics. An accessory kit contains five interchangeable chestpiece fittings – an adult and pediatric diaphragm, and three bells (1¼", 1", ¾") and deluxe eartips. A chestpiece key is also included to facilitate diaphragm replacement of the sealed gasket system. 22" tubing length. Five year warranty. Latex-free.

P542202	Sprague Stethoscope	Black
P542204	Sprague Stethoscope	Navy
P542206	Sprague Stethoscope	Neon Orange
P542208	Sprague Stethoscope	Neon Pink
P542210	Sprague Stethoscope	Red
P542212	Sprague Stethoscope	Royal Blue

Professional

Stethoscopes



Nurse Stethoscope

The Pro Advantage nurse stethoscope features a diaphragm only nurse style chestpiece with a non-chill rim, lightweight adjustable binaurals, flexible PVC 22" Y tubing, and choice of eartips. One year warranty. Latex-free.

P542002	Nurse Stethoscope	Black
P542007	Nurse Stethoscope	Gray
P542012	Nurse Stethoscope	Navy
P542017	Nurse Stethoscope	Red
P542022	Nurse Stethoscope	Royal Blue
P542027	Nurse Stethoscope	Teal



Dual Head Stethoscope

The Pro Advantage dual head stethoscope features a combination diaphragm/bell chestpiece with a non-chill rim, lightweight adjustable binaurals, flexible PVC 22" Y tubing, and choice of eartips. One year warranty. Latex-free.

P542010	Dual Head Stethoscope	Black
P542015	Dual Head Stethoscope	Gray
P542020	Dual Head Stethoscope	Navy
P542025	Dual Head Stethoscope	Red
P542030	Dual Head Stethoscope	Royal Blue
P542035	Dual Head Stethoscope	Teal



Bowles Stethoscope

The Pro Advantage bowles stethoscope features a diaphragm only chestpiece with a traditional raised stem design, lightweight adjustable binaurals, flexible PVC 22" Y tubing, and choice of eartips. One year warranty. Latex-free.

P542029	Bowles Stethoscope	Black
---------	--------------------	-------



Pediatric Stethoscope

The Pro Advantage pediatric stethoscope features a combination chestpiece with 1 1/4" diaphragm and 3/4" non-chill bell and rim, lightweight adjustable binaurals, flexible PVC 22" Y tubing, and choice of eartips. One year warranty. Latex-free.

P542111	Pediatric Stethoscope	Black
P542121	Pediatric Stethoscope	Gray



Infant Stethoscope

The Pro Advantage infant stethoscope features a combination chestpiece with a 25mm diaphragm and bell, slip on gray non-chill bell ring, lightweight adjustable binaurals, flexible PVC 22" Y tubing, and choice of eartips. One year warranty. Latex-free.

P542131	Infant Stethoscope	Black
P542141	Infant Stethoscope	Gray

Distributed By:

