



Analog & Digital TENS Unit

- Transcutaneous electrical nerve stimulation (TENS) therapy is the most common form of electrical stimulation used for pain management
- Small, battery-operated device delivers low-voltage electrical current through the skin via Pro Advantage electrodes placed near the source of pain
- Electricity from the electrodes stimulates nerves in the affected area and sends signals to the brain that “scramble” normal pain perception
- Available in analog and digital options



Item #	Description
P642000	2-Channel 3 Mode Digital TENS Unit
P641000	2-Channel Analog TENS Unit with Timer

Distributed by:

Gentle Stim Select

Neurostimulation Electrodes

- Clinically designed for all neurostimulation modalities such as TENS and muscle stimulation
- Strong tacky adhesive hydro-gel to meet the needs of patients
- High quality film for uniform dispersion over the electrode surface for maximum stimulation and comfort
- Low impedance
- Latex Free
- Pigtail end is pliable to fit most machine lead wires and designed for multiple uses
- Available in multiple shapes and sizes



Square



Round



Rectangle

Item #	Description	UOM
P640851	Cloth, 2" x 2" Square	4/pk, 10pk/bg
P640852	Cloth, 2" Round	4/pk, 10pk/bg
P640853	Cloth, 3" Round	4/pk, 10pk/bg
P640854	Cloth, 2" x 3.5" Rectangle	4/pk, 10pk/bg
P640855	Cloth, 1.5" x 3.5" Rectangle	4/pk, 10pk/bg
P640951	Foam, 2" x 2" Square	4/pk, 10pk/bg
P640952	Foam, 2" Round	4/pk, 10pk/bg
P640953	Foam, 3" Round	4/pk, 10pk/bg
P640954	Foam, 2" x 3.5" Rectangle	4/pk, 10pk/bg
P640955	Foam, 1.5" x 3.5" Rectangle	4/pk, 10pk/bg

Distributed by: