

Safety

Lancets

Used to safely obtain a capillary blood sample. Available in three color-coded versions to meet the desired blood sample size and address a variety of skin types.

- Single-use device
- Easy-to-use design ensures adequate blood flow at minimum pain level
- Needle is concealed before and after use to prevent needle-stick injuries and cross contamination exposure
- Ultra-sharp needle with low triggering force ensures gentle contact with puncture site
- Built-in needle depth control provides for consistent and accurate sampling results



Recommended use:

- 21G:**  **Normal Flow** For obtaining larger blood samples for testing glucose, hemoglobin, cholesterol, blood type studies, coagulation, blood gas and others
- 25G:**  **Low Flow** For blood glucose monitoring purposes. Can be used with older glucose monitors that require more blood.
- 28G:**  **Micro Flow** For glucose testing on children or patients with delicate skin. Ideal for modern glucometers that only require between 0.3ml – 1ml of blood.

Item #	Needle	Depth	Color	UOM
P929021	21G	1.8 mm	Blue	200/box
P929025	25G	1.5 mm	Purple	200/box
P929028	28G	1.2 mm	Green	200/box

Distributed By: