


DigiDop II



XD **eXtended** **Depth Probes**

56% more sensitive where you need it
- like shallow early fetal hearts! -
US PATENT PENDING

BREAKTHROUGH **SENSITIVITY** **XD – EXTENDED DEPTH PROBES**

www.Newman-Medical.com

AN HONEST

5 YEAR WARRANTY

Our 5 year warranty covers ALL labor and ALL parts Including probes and cords!



RIGHT SIZED

Extensive testing and feedback ensures the DigiDop fits comfortably in YOUR hand!



ASK YOUR SALES REPRESENTATIVE TODAY!
FOR QUESTIONS OR SALES SUPPORT CALL

800-267-5549

Smart Recharge System

Our premium rechargeable batteries stay charged for years when not in use and can be recharged more than 2000 times. A *significant* cost savings and reduction in waste.



DD-990R

Table top or wall mount version



PROBES	DESCRIPTION
D3	General purpose 3MHz probe - best for early fetal hearts
D2	2MHz probe designed for deeper fetal hearts
D3W	3MHz obstetrical for water applications
D2W	2MHz obstetrical for water applications

TECHNICAL SPECIFICATIONS	
Weight (with probe)	9.2 oz. (260 grams)
Dimensions	4.9x2.3x1.5in (125x60x40mm)
Warranty	5 years - all parts & labor
Materials	ABS/PC, stainless steel
Charging	Smart Recharge System
Probes	Interchangeable Probes
Charging Interface	Micro USB at 5V
Speaker Output	f0 260 Hz, 4W max
Country of Origin	United States of America

MAIN UNITS	DESCRIPTION
DD-990R	DigiDop II 990R, Table Top Doppler
DD-770AR	DigiDop II 770AR, Display w recharger & audio recorder
DD-770R	DigiDop II 770R, Display w recharger
DD-770	DigiDop II 770, Display
DD-330AR	DigiDop II 330AR, non-display w recharger & audio recorder
DD-330R	DigiDop II 330R, non-display w recharger
DD-330	DigiDop II 330, non-display



HOW TO ORDER

For product ordering code, add the probe number to the end of the main unit number with a hyphen. For example, to order a Display Doppler (DD-770) with 3MHz probe (D3) the order code is DD-770-D3.