

# Transport Chairs & Wheelchairs

2011  
nova



as unique as you are



Nova Transport Chairs and Wheelchairs bring ***Life in Motion*** with the best features and quality in the mobility world.

## FRAME STRENGTH

Nova Transport Chairs and Wheelchairs are made with the highest quality materials and workmanship. They are backed by our Limited Lifetime Warranty.

## BEST BRAKES

Whether it is the locking hand brakes on some models or wheel locks on others - Nova has the best brakes around!

## UPHOLSTERY

You can see and definitely feel the difference with Nova's seating upholstery. The material is attractive and comfortable, the craftsmanship is excellent, and the quality will last through life's many adventures.

## FOOTRESTS

All Nova Transport Chairs and Wheelchairs come with swing away, removable footrests. They are easy to adjust and durable for resting the unique shapes and sizes of feet out there.

## THE EXTRAS

With a Nova, you get the Extras as standard that come at an additional cost with other brands. Our 7000 Lightweight Wheelchairs come standard with Anti-Tippers and all Transport Chairs come standard with seat belts.

## **Tips for choosing your Nova Transport Chair or Wheelchair**

### **What's the main difference between a Wheelchair & a Transport Chair?**

Wheelchairs can be self-propelled by the user and also pushed by the companion. The large 24" rear wheels with rims allow the user to self-propel. This feature does add more weight and bulk to the wheelchair. Transport chairs have smaller rear wheels (8" & 12") so the user cannot self-propel and can only be pushed by the companion. Thus, the transport chair is both lighter and more compact than a wheelchair.

### **The Brakes - Who's in Control?**

#### **(Option for Transport Chairs)**

Some Nova Transport Chairs come with both rear wheel locks AND locking hand brakes at the companion handles. This allows the companion to control both the speed of the ride and ability to park the brake. This is a great option (& promotes safety) for those bumpy or hilly terrains. Make sure to have your seat belt on!

### **Weight and Width Matter!**

Be sure to understand and match the weight capacity and the seat width of the transport chair or wheelchair to the weight and general size of the user. You should also carefully consider the actual weight of the wheelchair or transport chair, especially if there are lots of events and adventures in the plans.

### **RIDE IN STYLE**

Choose the color and options that compliment your own personal tastes and lifestyle. Add a leopard Mobility Clutch for your cell phone, a flashlight attachment for the evening, a butterfly saddle bag for your important belongings, and a cup holder for your favorite beverages!

# STEEL TRANSPORT CHAIRS

## 307 & 309



### COMES STANDARD

- Fixed, full-length arms
- Carbon Steel Frame
- Quick release fold-down back
- Padded upholstery
- Padded armrests
- Locking rear wheels
- Removable aluminum footrests
- Non-skid foot plates & heel loops
- Foldable
- Includes seat belt
- 8" Wheels
- Seat width: 16.5"
- Weight: 26 lbs
- Weight Capacity: 300 lbs

### OPTIONS

- Travel Bag
- Cup Holder
- Flashlight Attachment
- Designer Mobility Bags
- Color: Hammertone





### COMES STANDARD

- Fixed, full-length arms
- Carbon Steel Frame
- Quick release fold-down back
- Padded upholstery
- Padded armrests
- Locking rear wheels
- Removable aluminum footrests
- Non-skid foot plates & heel loops
- Foldable
- Includes seat belt
- 8" Wheels
- Seat width: 18.5"
- Weight: 26 lbs
- Weight Capacities: 300 lbs

### OPTIONS

- Travel Bag
- Cup Holder
- Flashlight Attachment
- Designer Mobility Bags
- Colors: Black, Blue, Red



# ALUMINUM TRANSPORT CHAIRS

## 327 & 329



### COMES STANDARD

- Fixed, full-length arms
- Lightweight Aluminum frame
- Quick release fold-down back
- Padded upholstery
- Padded armrests
- Locking rear wheels
- Removable aluminum footrests
- Non-skid foot plates & heel loops
- Foldable
- Includes seat belt
- 8" Wheels
- Seat width: 16.5" (327) & 18.5" (329)
- Weight: 19 lbs
- Weight Capacities: 300 lbs

### OPTIONS

- Travel Bag
- Cup Holder
- Flashlight Attachment
- Designer Mobility Bags
- 327 Colors: Blue, Red
- 329 Colors: Blue, Red, Black





### COMES STANDARD

- Fixed, full-length arms
- Lightweight Aluminum Frame
- Quick release fold-down back
- Padded upholstery
- Padded armrests
- Locking rear wheels
- Removable aluminum footrests
- Non-skid foot plates & heel loops
- Foldable
- Includes seat belt
- 8" Wheels
- Seat width: 18"
- Weight: 18 lbs
- Weight Capacity: 250 lbs

### OPTIONS

- Travel Bag
- Cup Holder
- Flashlight Attachment
- Designer Mobility Bags
- Colors: Blue, Pink, Red



# ALUMINUM TRANSPORT CHAIRS

## 330



### COMES STANDARD

Lightweight Aluminum frame  
Hand Brakes for safety/control  
Quick release fold-down back  
Fixed, full-length arms  
Padded upholstery  
Padded armrests  
Locking rear wheels  
Removable aluminum footrests  
Non-skid foot plates & heel loops  
Includes seat belt  
Foldable  
12" Rear Wheels, 8" Front Wheels  
Seat width: 20.5"  
Weight: 25 lbs  
Weight Capacity: 300 lbs

### OPTIONS

Travel Bag  
Cup Holder  
Flashlight Attachment  
Designer Mobility Bags  
Colors: Blue, Red



***Dual Crossbars  
for Added Frame  
Strength!***

**332**



## COMES STANDARD

Lightweight Aluminum frame  
Hand Brakes for safety/control  
Detachable/Reversible desk  
arms  
Padded upholstery  
Padded armrests  
Locking rear wheels  
Removable aluminum footrests  
Non-skid foot plates & heel loops  
Includes seat belt  
12" Rear Wheels, 8" Front Wheels  
Seat width: 21.5"  
Weight: 27 lbs  
Weight Capacity: 400 lbs

## OPTIONS

Travel Bag  
Cup Holder  
Flashlight Attachment  
Designer Mobility Bags  
Colors: Blue, Red



# ALUMINUM TRANSPORT CHAIRS

349



## COMES STANDARD

Lightweight Aluminum frame  
Quick release fold-down back  
Detachable/Reversible desk  
arms  
Padded upholstery  
Padded armrests  
Locking rear wheels  
Removable aluminum footrests  
Non-skid foot plates & heel loops  
Foldable  
Includes seat belt  
8" Wheels  
Seat width: 18.5"  
Weight: 23 lbs  
Weight Capacity: 300 lbs

## OPTIONS

Travel Bag  
Cup Holder  
Flashlight Attachment  
Designer Mobility Bags  
Colors: Blue, Red





### COMES STANDARD

Lightweight Aluminum frame  
 Hand Brakes for safety/control  
 Quick release fold-down back  
 Detachable/Reversible desk arms  
 Padded upholstery  
 Padded armrests  
 Locking rear wheels  
 Removable aluminum footrests  
 Non-skid foot plates & heel loops  
 Includes seat belt  
 12" Rear Wheels, 8" Front Wheels  
 Seat width: 18.5"  
 Weight: 28 lbs  
 Weight Capacity: 300 lbs

### OPTIONS

Travel Bag  
 Cup Holder  
 Flashlight Attachment  
 Designer Mobility Bags  
 Colors: Blue, Red



# ALUMINUM TRANSPORT CHAIRS

## 378-R



### COMES STANDARD

Lightweight Aluminum frame  
Quick release fold-down back  
Fixed desk arms  
Padded armrests  
Locking rear wheels  
Removable, swing-away  
aluminum footrests  
Non-skid foot plates & heel loops  
Foldable  
Includes seat belt  
Removable 8" Wheels  
Travel bag included  
Seat width: 18"  
Weight: 19 lbs  
Weight Capacity: 250 lbs

### OPTIONS

Cup Holder  
Flashlight Attachment  
Designer Mobility Bags  
Colors: Blue, Red



379



## COMES STANDARD

Lightweight Aluminum frame  
Fixed desk arms  
Quick release fold-down back  
Locking rear wheels  
Removable, swing-away  
aluminum footrests  
Padded armrests  
Includes seat belt  
8" Wheels  
Seat width: 18.5"  
Weight: 28 lbs  
Weight Capacity: 300 lbs

## OPTIONS

Travel Bag  
Cup Holder  
Flashlight Attachment  
Designer Mobility Bags  
Colors: Blue, Red



# ALUMINUM TRANSPORT CHAIRS

368



## COMES STANDARD

- Lightweight aluminum frame
- Quick release drop arms
- Includes carry bag & strap
- Padded armrests
- Locking rear wheels
- Removable aluminum footrests
- Non-skid foot plates & heel loops
- Foldable
- Includes seat belt
- 8" Wheels
- Seat width: 18"
- Weight: 20.5 lbs
- Weight Capacity: 225 lbs

## OPTIONS

- Travel Bag
- Cup Holder
- Flashlight Attachment
- Designer Mobility Bags
- Colors: Blue, Red



TRANSPORT CHAIR SPECIFICATIONS

SPECIFICATIONS

Frame Model	307	309/319	327	329	348	330	332	349	352	368	378-R	379
Weight Capacity (lbs.)	300	300	300	300	250	300	400	300	300	225	250	300
Weight (lbs.)	26	26	19	19	18	25	27	23	28	20.5	19	17
Seat Width	16.5"	18.5"	16.5"	18.5"	18"	20.5"	21.5	18.5"	18.5"	18"	18"	18"
Seat Height from Floor	19.5"	19.5"	19.5"	19.5"	19"	19.5"	19.5"	19.5"	19.5"	19"	19"	19"
Arm Height from Seat	10"	10"	10"	10"	8"	12"	10"	10"	10"	10"	9"	9"
Back Height from Seat	20"	20"	20"	20"	17.5"	19"	17.5"	18"	18"	18"	17"	17"
Overall Height - Open	40.5"	40.5"	40.5"	40.5"	37"	40.5"	38"	38.25"	38"	29.25"	37.5"	37.5"
Overall Width - Open	20.5"	22.5"	20.5"	22.5"	22.5"	23.5"	28.5"	25"	25.5"	22.5"	23"	23"
Overall Height - Folded	30"	30"	30"	30"	28.5"	32"	38"	30.5"	30.5"	29.25"	30"	30"
Overall Width - Folded	11"	11"	11"	11"	12"	11.5"	11.5"	12"	13.5"	13"	9"	9"
Overall Depth - Folded	26"	26"	26"	26"	28"	27.5"	28.5"	26.5"	28"	13"	24"	24"

# TRANSPORT CHAIR ACCESSORIES



## CUP HOLDERS

- Item # CH-1000
- Item # CH-2000
- Attach easily to all walkers, wheelchairs, bed rails and transport chairs
- Adjust to hold different sized items (CH-1000 only)
- Folds up when not in use (CH-1000 only)



## TRAVEL BAG

- Item # 4000TB
- Nylon bag fits all Nova



## FLASHLIGHT

- Item # FL-1000
- Attaches easily to all walkers, wheelchairs, bed rails and transport chairs



## MOBILITY WRISTLETS

- Item # BGM01 (Black)
- Item # BGM04 (Zebra)
- Item # BGM03 (Pink)
- Item # BGM05 (Metropolitan)



## MOBILITY CLUTCHES

- Item # BGS01 (Black)
- Item # BGS02 (Blue)
- Item # BGS03 (Red)
- Item # BGS04 (Leopard)
- Item # BGS05 (Camouflage)

- Attach securely to all walkers, wheelchairs, transport chairs, and bed rails

## MOBILITY SADDLE BAGS



- Attach securely to all walkers, wheelchairs, transport chairs, and bed rails



- Item # 4002AP (Aloha Pink)
- Item # 4002BF (Butterflies)
- Item # 4002BK (Black)
- Item # 4002CT (Cats)
- Item # 4002DG (Dogs)
- Item # 4002MP (Metropolitan)

## 5060S/5080S

## K3



### 16"/18" WIDTHS

5060S/5080S - Full arms with s/a footrests

5060SE/5080SE - Full arms with elevating leg rest

### FEATURES

Dual axle positions for quick seat height conversion

3-position front castor forks

High strength aluminum foot plates

"Easy Adjust" foot riggings

Sealed wheel bearings

Flame retardant black vinyl upholstery

Industry compatible accessories

Back pocket

Lifetime Limited Warranty on Frame

Weights: 37 lbs

Weight Capacities: 275 lbs



# HAMMERTONE WHEELCHAIRS

## 5060H/5080H

## K3



### 16"/18"WIDTHS

5060H - Fixed arms and s/a footrests

5080H - Fixed arms and s/a footrests

### FEATURES

Dual axle positions for quick seat height conversion

3-position front castor forks

High strength aluminum foot plates

"Easy Adjust" foot riggings

Sealed wheel bearings

Flame retardant black vinyl upholstery

Industry compatible accessories

Back pocket

Lifetime Limited Warranty on Frame

Weights: 37 lbs

Weight Capacities: 275 lbs



## 5090H

## K1



### 20" SEAT WIDTH

5090H - Fixed arms and s/a footrests

### FEATURES

Dual axle positions for quick seat height conversion

3-position front castor forks

High strength aluminum foot plates

"Easy Adjust" foot riggings

Sealed wheel bearings

Flame retardant black vinyl upholstery

Industry compatible accessories

Back pocket

Lifetime Limited Warranty on Frame

Weight: 37 lbs

Weight Capacity: 275 lbs



# STEEL WHEELCHAIRS

## 5165S/5185S

## K1



### 16"/18" WIDTHS

**5165S/5185S** - Desk arms with s/a footrests

**5165SE/5185SE** - Desk arms with elevating leg rest

**5166S/5186S** - Full arms with s/a footrest

**5166SE/5186SE** - Full arms with elevating leg rest

### FEATURES

Detachable/Reversible desk arms

Dual axle positions for quick seat height conversion

3-position front castor forks

High strength aluminum foot plates

"Easy Adjust" foot riggings

Sealed wheel bearings

Flame retardant black vinyl upholstery

Industry compatible accessories

Lifetime Limited Warranty on Frame

Weights: 42 lbs

Weight Capacities: 275 lbs



## 5160S/5180S

## K3



### 16"/18" WIDTHS

**5160S/5180S** - Desk arms with s/a footrests

**5160SE/5180SE** -Desk arms with elevating leg rest

**5161S/5181S** - Full arms with s/a footrest

**5161SE/5181SE** - Full arms with elevating leg rest

### FEATURES

Detachable/Reversible desk arms

Dual axle positions for quick seat height conversion

3-position front castor forks

High strength aluminum foot plates

"Easy Adjust" foot riggings

Sealed wheel bearings

Flame retardant black nylon upholstery

Industry compatible accessories

Lifetime Limited Warranty on Frame

Weights: 35 lbs

Weight Capacities: 275 lbs



# STEEL WHEELCHAIRS

## 5200S

## K3



### 20" WIDTH

**5200S** - Desk arms with s/a footrests

**5200SE** - Desk arms with elevating leg rest

**5201S** - Full arms with s/a footrest

**5201SE** - Full arms with elevating leg rest

### FEATURES

Detachable/Reversible desk arms

Dual axle positions for quick seat height conversion

3-position front castor forks

High strength aluminum foot plates

"Easy Adjust" foot riggings

Sealed wheel bearings

Flame retardant black nylon upholstery

Industry compatible accessories

Lifetime Limited Warranty on Frame

Weight: 35 lbs

Weight Capacity: 325 lbs



## 5220S/5240S

## K7



***Dual Crossbars  
for Added Frame  
Strength!***

### 22"/24" WIDTHS (18" SEAT DEPTH)

5220S/5240S - Desk arms with s/a footrests

5220SE/5240SE - Desk arms with elevating leg rest

5221S/5241S - Full arms with s/a footrest

5221SE/5241SE - Full arms with elevating leg rest

### FEATURES

Dual axle positions for quick seat height conversion

3-position front castor forks & dual cross bar supports

High strength aluminum foot plates

"Easy Adjust" foot riggings

Sealed wheel bearings

Flame retardant black vinyl upholstery

Industry compatible accessories

Back pocket

Lifetime Limited Warranty on Frame

Weights (5220/5221): 56 lbs

Weights (5240/5241): 57 lbs

Weight Capacities (5220/5221): 350 lbs

Weight Capacities (5240/5241): 425 lbs



# RECLINING WHEELCHAIRS

## 6160S/6180S/6200S

## K1



### 16"/18"/20" WIDTHS

**6160S** - Recliner, detachable full arms, elevating leg rest

**6180S** - Recliner, detachable full arms, elevating leg rest

**6200S** - Recliner, detachable full arms, elevating leg rest

### FEATURES

Reclines to full 180 degree position

"Easy Adjust" foot riggings

Sealed wheel bearings

High strength aluminum foot plates

Flame retardant black vinyl upholstery

Industry compatible accessories

Lifetime Limited Warranty on Frame

Weight (6160S): 46 lbs

Weights (6180S/6200S): 48 lbs

Weight Capacities (6160S/6180S): 275 lbs

Weight Capacity (6200S): 325 lbs



## 7160L/7180L

## K4



### 16"/18" WIDTHS

7160L/7180L - Flip back desk arms with s/a footrests

7160LE/7180LE - Flip back desk arms w/ elevating leg rests

7161L/7181L - Full arms with s/a footrests

7161LE/7181LE - Flip back desk arms w/ elevating leg rests

### FEATURES

Dual axle positions for quick seat height conversion

Flip back desk arms

3-position front castor forks

High strength aluminum foot plates

"Easy Adjust" foot riggings

Flame retardant black nylon upholstery

Includes heel loops and anti-tippers

Industry compatible accessories

Lifetime Limited Warranty on Frame

Weights: 32 lbs

Weight Capacities: 300 lbs



# LIGHTWEIGHT WHEELCHAIRS

## 7190L

## K4



### 18" WIDTH (18" SEAT DEPTH)

7190L - Flip back desk arms with s/a footrests

7190LE - Flip back desk arms w/ elevating leg rests

### FEATURES

Dual axle positions for quick seat height conversion

Dual crossbar supports

Flip back desk arms

3-position front castor forks

High strength aluminum foot plates

"Easy Adjust" foot riggings

Flame retardant black nylon upholstery

Industry compatible accessories

Lifetime Limited Warranty on Frame

Weight: 46 lbs

Weight Capacity: 300 lbs



## 7200L

## K4



### 20" WIDTH

**7200L** - Flip back desk arms with s/a footrests

**7200LE** - Flip back desk arms w/ elevating leg rests

**7201L** - Full arms with s/a footrests

**7201LE** - Flip back desk arms w/ elevating leg rests

### FEATURES

Dual axle positions for quick seat height conversion

Dual crossbar supports

Flip back desk arms

3-position front castor forks

High strength aluminum foot plates

"Easy Adjust" foot riggings

Flame retardant black nylon upholstery

Industry compatible accessories

Lifetime Limited Warranty on Frame

Weight : 33 lbs

Weight Capacity: 300 lbs



# WHEELCHAIR SPECIFICATIONS

## SPECIFICATIONS

*\* Includes Headrest*

Frame Model	5060S	5080S	5060H	5080H	5090H	5165S	
Weight Capacity (lbs.)	275	275	275	275	275	275	
Weight (lbs.)	37	37	37	37	37	42	
Seat Width	16"	18"	16"	18"	20"	16"	
Seat Height from Floor	17.5"-19.5"	17.5"-19.5"	17.5"-19.5"	17.5"-19.5"	17.5"-19.5"	17.5"-19.5"	
Arm Height from Seat	10"	10"	10"	10"	10"	8.5"	
Back Height from Seat	15"	15"	15"	15"	15"	15"	
Overall Height	35.5"	35.5"	35.5"	35.5"	35.5"	35.5"	
Overall Width	22.5"	24.5"	22.5"	24.5"	26.5"	23.5"	
Overall Depth	31.5"	31.5"	31.5"	31.5"	31.5"	31.5"	
Overall Width - Folded	11"	11"	11"	11"	11"	12"	
Seat Depth	16"	16"	16"	16"	16"	16"	
Frame Model	5185S	5160S	5180S	5200S	5220S	5240S	
Weight Capacity (lbs.)	275	275	275	325	350	425	
Weight (lbs.)	42	35	35	35	56	57	
Seat Width	18"	16"	18"	20"	22"	24"	
Seat Height from Floor	17.5"-19.5"	17.5"-19.5"	17.5"-19.5"	17.5"-19.5"	17.5"-19.5"	17.5"-19.5"	
Arm Height from Seat	8.5"	9"	9"	9"	9"	9"	
Back Height from Seat	15"	15"	15"	16"	16"	16"	
Overall Height	35.5"	35.5"	35.5"	35.5"	36.25"	36.25"	
Overall Width	25.5"	24"	26"	29"	31"	33"	
Overall Depth	31.5"	31.5"	31.5"	31.5"	34"	35"	
Overall Width - Folded	12"	12"	12"	12"	13.5"	13.5"	
Seat Depth	16"	16"	16"	16"	18"	18"	
Frame Model	6160S	6180S	6200S	7160L	7180L	7190L	7200L
Weight Capacity (lbs.)	275	275	325	300	300	300	300
Weight (lbs.)	46/*49	48/*51	48/*51	32	32	46	33
Seat Width	16"	18"	20"	16"	18"	18"	20"
Seat Height from Floor	20"	20"	20"	17.5"-19.5"	17.5"-19.5"	17.5"-19.5"	17.5"-19.5"
Arm Height from Seat	8.5"	8.5"	8.5"	10"	10"	10"	10"
Back Height from Seat	16"	16"	16"	16"	16"	16"	16"
Overall Height	37"/*49"	37"/*49"	37"/*49"	36.25"	36.25"	36.25"	36.25"
Overall Width	24.25"	26.25"	28.25"	23"	25"	25"	27"
Overall Depth	36.5"	36.5"	36.5"	31.5"	31.5"	31.5"	31.5"
Overall Width - Folded	12.5"	12.5"	12.5"	11"	11"	11"	11"
Seat Depth	16"	16"	16"	16"	16"	18"	16"

**NEW!**



## Glamour Line



### ULTIMATE WALKER BAGS

• Versatile to fit on rolling walkers and folding walkers. Generous and roomy with convenient front pockets.

- Item # 4003BK (Black)
- Item # 4003EA (Edgy Aqua)
- Item # 4003HT (Houndstooth)
- Item # 4003LD (Leopard)

### UNIVERSAL MOBILITY BAGS

• This bag goes global! Fits on any bar, backrest or handles. So many options, so many places to go. Large main compartment, multiple outer pockets.

- Item # 4004BK (Black)
- Item # 4004EA (Edgy Aqua)
- Item # 4004HT (Houndstooth)
- Item # 4004LD (Leopard)

### MOBILITY HAND BAGS

• So sleek and stylish. Perfect to go from walker or rail to over the shoulder. (Comes with adjustable shoulder strap).

- Item # 4005BK (Black)
- Item # 4005EA (Edgy Aqua)
- Item # 4005HT (Houndstooth)
- Item # 4005LD (Leopard)

### SEAT AND BACK COVERS

• Add glam, add style, add luxurious feel. Make it all about you. Stand out from the crowd and make a statement!

- Item # 4007SC (Safari Cheetah)
- Item # 4007SL (Snow Leopard)

# WHEELCHAIR CUSHIONS

## GEL FOAM WHEELCHAIR CUSHION

<b>ITEM # 2601-2-1616</b>	2" Gel Foam Cushion for 16" x 16" Wheelchair
<b>ITEM # 2601-2-1816</b>	2" Gel Foam Cushion for 18" x 16" Wheelchair
<b>ITEM # 2601-3-1616</b>	3" Gel Foam Cushion for 16" x 16" Wheelchair
<b>ITEM # 2601-3-1816</b>	3" Gel Foam Cushion for 18" x 16" Wheelchair
<b>ITEM # 2601-3-2016</b>	3" Gel Foam Cushion for 20" x 16" Wheelchair


- Gel/foam combination provides pressure relief, prevents minor skin irritations, and improves weight distribution
- Removable, machine washable polyester/cotton cover
- Non-skid bottom with velcro straps secure cushion to wheelchair



## CONVOLUTED SEAT/BACK CUSHION

<b>ITEM # 2657-3</b>	3" Convoluted Seat/Back Foam Cushion without Cover for 18" x 16" Wheelchair - Blue
<b>ITEM # 2658-3</b>	3" Convoluted Seat/Back Foam Cushion with Fleece Cover for 18" x 16" Wheelchair - Blue

- Provides support and pressure relief
- Convoluted foam improves weight distribution and air circulation
- 250 lb. weight capacity
- 1 year limited warranty against defects in workmanship
- Removable, machine washable polyester/cotton for added comfort (Fleece Cover)



# WHEELCHAIR CUSHIONS

## 3" COCCYX FOAM CUSHION

<b>ITEM # 2603C-3-1616</b>	For 16" x 16" Wheelchair
<b>ITEM # 2603C-3-1816</b>	For 18" x 16" Wheelchair
<b>ITEM # 2603C-3-2016</b>	For 20" x 16" Wheelchair
<b>ITEM # 2603C-3-2018</b>	For 20" x 18" Wheelchair


- Coccyx cutout reduces pressure on tailbone and spine
- Cushion helps equalize pressure and reduce soreness
- Waterproof cover is removable and machine washable
- 250 lb. weight capacity
- 1 year limited warranty against defects in workmanship



## 3" CONVOLUTED FOAM CUSHION

<b>ITEM # 2650-3</b>	Without Cover for 18" x 16" Wheelchair
<b>ITEM # 2651-3</b>	With Cover for 18" x 16" Wheelchair
<b>ITEM # 2652-3</b>	With Cover for 16" x 16" Wheelchair

- Provides support and pressure relief
- Slip free seat promotes proper positioning and improves posture
- Convoluted foam improves weight distribution and air circulation
- Highly resilient convoluted polyurethane foam provides comfort and durability



# WHEELCHAIR CUSHIONS

## Back Foam Wheelchair Cushion with Lumbar Support and Stabilization Board

ITEM # 2611BK-16	16" Back foam cushion
ITEM # 2611BK-18	18" Back foam cushion
ITEM # 2611BK-20	20" Back foam cushion

- Curved profile fits natural back contours
- Non-flammable and cover is washable and waterproof
- Stabilization board prevents wheelchair sling back from "hammocking"
- Nylon straps securely position cushion to wheelchair
- Weight Capacity: 250 lbs
- 1 Year Limited Warranty against defects in workmanship
- Dimensions: 2.5"h x 16"w x 16"d (for 16" wide)  
2.5"h x 18"w x 16"d (for 18" wide)



# WHEELCHAIR CUSHIONS

## Foam Coccyx Cushion

**ITEM # 2654C-3** For 18" x 16" Wheelchair

- Coccyx cutout reduces pressure on the tailbone and spine
- Wheelchair cushion helps equalize pressure and reduce soreness
- Removable, machine washable polyester/cotton cover
- Made of High-Density Foam
- Flame retardant material (CA 117)
- Weight Capacity: 250 lbs
- 1 Year Limited Warranty against defects in workmanship



## 3" Foam Wheelchair Cushion

**ITEM # 2664-3** With Cover for 18" x 16" Wheelchair

- Wheelchair cushion helps equalize pressure and reduce soreness
- Removable, machine washable polyester/cotton cover
- Made of High-Density Foam
- Fits 18" x 16" Wheelchair
- Flame retardant material (CA 117)
- Weight Capacity: 250 lbs
- 12 Month Limited Warranty against defects in workmanship



# MORE ACCESSORIES!

## ANTI-TIPPERS

- Helps prevent wheelchairs from tipping backward.



**Item # Fits/Serial Numbers/Color**

- **W-5105** • 5060S,5080S,5165S,5185S,5220S,5240S/Y,J,H/Chrome
- **W-5106** • 5060S,5080S,5165S,5185S,5220S,5240S/Y,J,H/Chrome
- **W-5005** • 5200S,5180S,7160S,7200L,7180L,7160L/Y,J,H/Black
- **J-AT100** • 5200S,5180S,7160S,7200,7180L,7160L/Y,J,H/Chrome
- **J-AT101** • 5060S,5080S,5165S,5185S,5220S,5240S/Y,J,H/Chrome

## SEAT BELT

- **Item # W-7006**
- Airline Style



## BRAKE EXTENSIONS

- **Item # J-BE100**
- For 7000 Series Wheelchairs

## HANGING MOBILITY POUCH

- Convenient for carrying various personal belongings
- Easily attaches to all Wheelchairs, Transport Chairs, Walkers and Bed Rails



• **Item # 4001WP**

• **Item # P4200SB**



# Service



Eighteen years ago, we made the commitment to put our customers first. We honor our warranties and we care. Simply put, we will do “whatever it takes.”

NOVA is a leading manufacturer of healthcare Mobility and Bathroom Safety products. Internationally started in 1981 as a family owned and operated business, our vision was to provide quality home care products as a means of promoting a healthy and independent lifestyle. NOVA remains family owned and operated and has earned a worldwide reputation for quality, integrity, and service.

NOVA’s commitment to customers begins by manufacturing each product with the highest quality standards. Our commitment continues for the life of the product with caring and professional customer support. All NOVA Transport Chairs and Wheelchairs are backed by a Limited Lifetime Warranty.

***We welcome you to NOVA’s  
Family of Healthcare Products!***



*Experience the NOVA Way!*  
[www.NovaMedicalProducts.com](http://www.NovaMedicalProducts.com)

**LIVE LIFE TO THE  
FULLEST WITH YOUR NOVA:  
The very best in Transport Chairs &  
Wheelchairs.**



Improve your overall mobility, safety, independence, and style by selecting from Nova's extensive line of Transport Chairs and Wheelchairs.

**WARNING**

- Lock wheel locks when sitting stationary
- Do not use on stairs, escalators, or moving walkways
- Consult your user manual before operating



**as unique as you are**

**1-800-557-6682**

**[www.NovaMedicalProducts.com](http://www.NovaMedicalProducts.com)**

Authorized NOVA Dealer