

Rolling Walkers



as unique as you are



Nova Rolling Walkers are **Built for Life** with the best features and quality in the mobility world.

FRAME STRENGTH

Nova rolling walker frames are made with the highest quality materials and workmanship. They are backed by our Limited Lifetime Warranty.

HAND BRAKES

Our **patented hand brakes** are designed for maximum comfort and control. They allow users to easily control the speed by pulling up and engage the park brake by pushing down. Users with weak or small hands need not worry – our hand brakes are the best and easiest to use. They even come with reflectors for safety and good looks!

FLIP-UP BACK

(On 4200, 4200C, 4201, 4202, 4202C, 4203, 4204, 4207 4208, 4212, 4214, 4215, 4218, 4220, 4222, 4224)

Nova's unique **Flip-up Back** is so convenient for compact transporting and just as easily flips back down to provide a comfortable backrest when the user is seated. No more lost or damaged backrests with Nova's easy Flip-up Back.

EASY TO REMOVE BACK

(On 4216)

Nova's **patented pull pin release design** allows the user to easily remove the back for compact storage. The back can just as easily lock back into place for a secure and comfortable backrest.

Tips for Choosing Your Nova Rolling Walker

Walkers are NOT one size fits all

To determine the correct seat height – sit on the seat with your back to the backrest. Your feet should both be flat on the ground with your thighs resting comfortably on the seat. Make sure your weight is below the listed weight capacity and you can sit comfortably between the handles.

Choose the features and options that go with your lifestyle

With your Nova walker you will now be on the go! - so choose the features that compliment and support your independence. What kind of seat do you want? How about storage for your personal belongings? What size wheels are best for your walking terrain? Be sure to consider great options like a cup holder, flashlight, or mobility clutch!

Get your Style ON

Choose the walker color and designer mobility bags to compliment your personal style. You can choose for YOU and walk with pride and style.

Before you take your Nova for a test drive – be sure to:

- Make sure the backrest is attached.
- Adjust the handles to your height. Adjust the hand grips to your wrist for a slight 30 degree bend at the elbow.
- Read and understand all user warnings.
- Practice engaging the hand brakes and locking the hand brakes prior to use.

Zoom 4218

Get the perfect fit and features with the **Zoom**. Choose your perfect seat height from **18," 20," 22,"** and **24"** for your height, from petite to tall. The **Zoom** has the best features with 8" wheels, a large padded seat, and an under seat pouch.



COMES STANDARD

- Locking Hand Brakes
- Large Padded Seat
- Under Seat Pouch
- Flip-up Back
- Foldable
- 8" Wheels

OPTIONS

- Travel Bag
- Cup Holder
- Flashlight Attachment
- Cane Holder
- Designer Mobility Bags
- Colors: Blue, Pink, Red

WEIGHT CAPACITY: 300 lbs

SEAT HEIGHT: 18"

WEIGHT: 15.5 lbs

SEAT DIMENSIONS: 13.5"w x 12.5"d

HANDLE HEIGHT ADJUSTMENT: 28" - 31.25"

WIDTH BETWEEN HANDLES: 17.5"

OVERALL DIMENSIONS: 24.5"w x 25"d

APPROX. USER HEIGHT: 4'10" - 5'4"

POUCH DIMENSIONS: 9"l x 15"w x 6.5"d

DIMENSIONS WHEN FOLDED: 34"l x 24.5"w x 10.5"d



Zoom 4220

Get the perfect fit and features with the **Zoom**. Choose your perfect seat height from **18," 20," 22,"** and **24"** for your height, from petite to tall. The **Zoom** has the best features with 8" wheels, a large padded seat, and an under seat pouch.



COMES STANDARD

- Locking Hand Brakes
- Large Padded Seat
- Under Seat Pouch
- Flip-up Back
- Foldable
- 8" Wheels

OPTIONS

- Travel Bag
- Cup Holder
- Flashlight Attachment
- Cane Holder
- Designer Mobility Bags
- Colors: Black, Blue, Green, Red

WEIGHT CAPACITY: 300 lbs

SEAT HEIGHT: 20"

WEIGHT: 16.5 lbs

SEAT DIMENSIONS: 13.5"w x 12.5"d

HANDLE HEIGHT ADJUSTMENT: 29.5" - 34"

WIDTH BETWEEN HANDLES: 17.5"

OVERALL DIMENSIONS: 24.5"w x 26"d

APPROX. USER HEIGHT: 5'2" - 5'8"

POUCH DIMENSIONS: 9"l x 15"w x 6.5"d

DIMENSIONS WHEN FOLDED: 34"l x 24.5"w x 10.5"d



Zoom 4222

Zoom 4224

Get the perfect fit and features with the **Zoom**. Choose your perfect seat height from **18," 20," 22,"** and **24"** for your height, from petite to tall. The **Zoom** has the best features with 8" wheels, a large padded seat, and an under seat pouch.

Get the perfect fit and features with the **Zoom**. Choose your perfect seat height from **18," 20," 22,"** and **24"** for your height, from petite to tall. The **Zoom** has the best features with 8" wheels, a large padded seat, and an under seat pouch.



COMES STANDARD

- Locking Hand Brakes
- Large Padded Seat
- Under Seat Pouch
- Flip-up Back
- Foldable
- 8" Wheels

OPTIONS

- Travel Bag
- Cup Holder
- Flashlight Attachment
- Cane Holder
- Designer Mobility Bags
- Colors: Black, Blue, Green, Red

COMES STANDARD

- Locking Hand Brakes
- Large Padded Seat
- Under Seat Pouch
- Flip-up Back
- Foldable
- 8" Wheels

OPTIONS

- Travel Bag
- Cup Holder
- Flashlight Attachment
- Cane Holder
- Designer Mobility Bags
- Colors: Blue, Red

WEIGHT CAPACITY: 300 lbs **WEIGHT:** 16.5 lbs
SEAT HEIGHT: 22" **SEAT DIMENSIONS:** 13.5"w x 12.5"d
HANDLE HEIGHT ADJUSTMENT: 31.5" - 36"
WIDTH BETWEEN HANDLES: 17.5"
OVERALL DIMENSIONS: 24.5"w x 26.5"d
APPROX. USER HEIGHT: 5'5" - 6'1"
BAG DIMENSIONS: 9"l x 15"w x 6.5"d
DIMENSIONS WHEN FOLDED: 34"l x 24.5"w x 10.5"d



WEIGHT CAPACITY: 300 lbs **WEIGHT:** 16.5 lbs
SEAT HEIGHT: 24" **SEAT DIMENSIONS:** 13.5"w x 12.5"d
HANDLE HEIGHT ADJUSTMENT: 33.5" - 37.75"
WIDTH BETWEEN HANDLES: 17.5"
OVERALL DIMENSIONS: 24"w x 27.5"d
APPROX. USER HEIGHT: 5'5" - 6'2"
BAG DIMENSIONS: 9"l x 15"w x 6.5"d
DIMENSIONS WHEN FOLDED: 36"l x 24.5"w x 10.5"d



GetGO Groove 4204

The **GetGO Groove** has all of the best rolling walker features, including **two seat heights**. It can easily be assembled and unassembled for convenient travel and storage.



COMES STANDARD

- Locking Hand Brakes
- Large Padded Seat
- Under Seat Pouch
- Flip-up Back
- Foldable/Collapsible
- 8" Wheels

OPTIONS

- Travel Bag
- Cup Holder
- Flashlight Attachment
- Cane Holder
- Designer Mobility Bags
- Colors: Blue, Red

WEIGHT CAPACITY: 250 lbs **WEIGHT:** 16.5 lbs

SEAT HEIGHT ADJUSTMENT: 19.5" - 21.5"

SEAT DIMENSIONS: 13.75"w x 12.75"d

HANDLE HEIGHT ADJUSTMENT: 31.75" - 35.25"

WIDTH BETWEEN HANDLES: 18"

OVERALL DIMENSIONS: 24.5"w x 27.5"d

APPROX. USER HEIGHT: 5'3" - 6'0"

BAG DIMENSIONS: 7"l x 15"w x 8.5"d

DIMENSIONS WHEN FOLDED: 42"l x 23"w x 8.5"d



GetGO Classic 4202C

The **GetGO Classic** is a true classic with **great features** at a **nice value**. Equipped with 6" wheels, a large padded seat and an under seat pouch.



COMES STANDARD

- Locking Hand Brakes
- Large Padded Seat
- Under Seat Pouch
- Flip-up Back
- Foldable
- 6" Wheels

OPTIONS

- Travel Bag
- Cup Holder
- Flashlight Attachment
- Cane Holder
- Designer Mobility Bags
- Colors: Blue, Red, Green, Pink

WEIGHT CAPACITY: 300 lbs

WEIGHT: 13.5 lbs

SEAT HEIGHT: 22"

SEAT DIMENSIONS: 13.5"w x 12.5"d

HANDLE HEIGHT ADJUSTMENT: 32.25" - 36.5"

WIDTH BETWEEN HANDLES: 17.5"

OVERALL DIMENSIONS: 24.5"w x 24"d

APPROX. USER HEIGHT: 5'4" - 6'1"

BAG DIMENSIONS: 9"l x 15"w x 6.5"d

DIMENSIONS WHEN FOLDED: 36"l x 24.5"w x 9.5"d



GetGO 4203

The **GetGO** will get you on the go with **great features** at a very **nice value**. Equipped with 6" wheels, a plastic seat with built-in seat handle, and a front wire basket.



COMES STANDARD

Locking Hand Brakes
Plastic Seat with built-in handle
Front Wire Basket
Flip-up Back
Foldable
6" Wheels

OPTIONS

Travel Bag
Basket Cover Bag
Tray
Rubber Seat Pad
Cup Holder
Flashlight Attachment
Cane Holder
Designer Mobility Bags
Colors: Black, Blue, Pink, Pink, Red

WEIGHT CAPACITY: 300 lbs

WEIGHT: 13.5 lbs

SEAT HEIGHT: 21.5"

SEAT DIMENSIONS: 13.5"w x 13"d

HANDLE HEIGHT ADJUSTMENT: 32.25" - 36.5"

WIDTH BETWEEN HANDLES: 17.5"

OVERALL DIMENSIONS: 24"w x 24"d

APPROX. USER HEIGHT: 5'4" - 6'1"

BASKET DIMENSIONS: 10"l x 15"w x 6.5"d

DIMENSIONS WHEN FOLDED: 36"l x 24.5"w x 8"d



GetGO Petite 4208

The **GetGO Petite**, ideal for petite and slight users will get you on the go with **great features** at a very **nice value**.



COMES STANDARD

Locking Hand Brakes
Plastic Seat with built-in handle
Front Wire Basket
Flip-up Back
Foldable
6" Wheels

OPTIONS

Travel Bag
Rubber Seat Pad
Cup Holder
Basket Cover Bag
Flashlight Attachment
Cane Holder
Designer Mobility Bags
Colors: Blue, Pink, Red

WEIGHT CAPACITY: 300 lbs

WEIGHT: 12 lbs

SEAT HEIGHT: 18"

SEAT DIMENSIONS: 13.5"w x 12.5"d

HANDLE HEIGHT ADJUSTMENT: 29" - 33"

WIDTH BETWEEN HANDLES: 15.5"

OVERALL DIMENSIONS: 22"w x 23.5"d

APPROX. USER HEIGHT: 4'10" - 5'4"

BASKET DIMENSIONS: 8.5"l x 13.5"w x 6.5"d

DIMENSIONS WHEN FOLDED: 33"l x 22"w x 7"d



The **GetGO Super** is ideal for taller users. It will get you on the go with **great features** at a **very nice value**.



COMES STANDARD

- Locking Hand Brakes
- Plastic Seat with built-in handle
- Front Wire Basket
- Flip-up Back
- Foldable
- 8" Wheels

OPTIONS

- Travel Bag
- Tray
- Rubber Seat Pad
- Cup Holder
- Basket Cover Bag
- Flashlight Attachment
- Cane Holder
- Designer Mobility Bags
- Colors: Blue, Red

WEIGHT CAPACITY: 300 lbs

WEIGHT: 16.5 lbs

SEAT HEIGHT: 23.5"

SEAT DIMENSIONS: 13.5"w x 13"d

HANDLE HEIGHT ADJUSTMENT: 34.5" - 39"

WIDTH BETWEEN HANDLES: 17.5"

OVERALL DIMENSIONS: 22"w x 28.5"d

APPROX. USER HEIGHT: 5'6" - 6'2"

BASKET DIMENSIONS: 9.5"l x 15"w x 6.5"d

DIMENSIONS WHEN FOLDED: 33"l x 22"w x 10"d



The **Traveler** will take you many places with its unique, **three wheel** maneuverable design.



COMES STANDARD

- Locking Hand Brakes
- Tray
- Basket
- Pouch
- Foldable
- 8" Wheels

OPTIONS

- Cup Holder
- Travel Bag
- Flashlight Attachment
- Cane Holder
- Designer Mobility Bags
- Colors: Black, Blue, Red

WEIGHT CAPACITY: 250 lbs

WEIGHT: 15 lbs

HANDLE HEIGHT ADJUSTMENT: 30.25" - 35.25"

WIDTH BETWEEN HANDLES: 12"

OVERALL DIMENSIONS: 25"w x 24"d

APPROX. USER HEIGHT: 5'4" - 6'0"

BASKET DIMENSIONS: 9"l x 10"w x 9.5"d

BAG DIMENSIONS: 13"l x 12"w x 9.5"d

TRAY DIMENSIONS: 9.25"w x 8"d

DIMENSIONS WHEN FOLDED: 27"l x 11"w x 26.5"d



Cruiser Deluxe 4202

The **Cruiser Deluxe** is the **legend** in the Nova Line. Equipped with the **best of the best**, the **Cruiser Deluxe** comes with 8" wheels, a deluxe seat with built-in seat pad, and a sleek frame.



COMES STANDARD

Locking Hand Brakes
Plastic seat with built-in rubber seat pad
Front Wire Basket
Removable Back
Foldable
8" Wheels

OPTIONS

Travel Bag
Tray
Basket Cover Bag
Seat Cover
Cup Holder
Flashlight Attachment
Cane Holder
Designer Mobility Bags
Colors: Black, Blue, Green, Purple, Red

WEIGHT CAPACITY: 300 lbs

WEIGHT: 17.5 lbs

SEAT HEIGHT: 21"

SEAT DIMENSIONS: 16.75"w x 12.5"d

HANDLE HEIGHT ADJUSTMENT: 32" - 36"

WIDTH BETWEEN HANDLES: 18.5"

OVERALL DIMENSIONS: 23"w x 28"d

APPROX. USER HEIGHT: 5'4" - 6'0"

BASKET DIMENSIONS: 8"l x 16.5"w x 6"d

DIMENSIONS WHEN FOLDED: 35"l x 23"w x 10"d



Cruiser Deluxe Junior 4207

The **Cruiser Deluxe Junior** has a lower handle adjustment and seat height than the Cruiser Deluxe. The **Cruiser Deluxe Junior** comes with 6" wheels, a deluxe seat with built-in pad, and a sleek frame.



COMES STANDARD

Locking Hand Brakes
Plastic seat with built-in rubber seat pad
Front Wire Basket
Removable Back
Foldable
6" Wheels

OPTIONS

Travel Bag
Tray
Basket Cover Bag
Seat Cover
Cup Holder
Flashlight Attachment
Cane Holder
Designer Mobility Bags
Colors: Blue, Purple, Red

WEIGHT CAPACITY: 300 lbs

WEIGHT: 16.5 lbs

SEAT HEIGHT: 19"

SEAT DIMENSIONS: 16.75"w x 12.5"d

HANDLE HEIGHT ADJUSTMENT: 30.25" - 34.5"

WIDTH BETWEEN HANDLES: 18.5"

OVERALL DIMENSIONS: 23"w x 26"d

APPROX. USER HEIGHT: 4'11" - 5'5"

BASKET DIMENSIONS: 8"l x 16"w x 6"d

DIMENSIONS WHEN FOLDED: 35"l x 23"w x 11"d



Cruiser 4200

The **Cruiser** allows you to cruise along with a **weight activated** brake system. This walker **DOES NOT** have locking hand brakes for those unable to use them. Moderate-to-good stability is required to use the non-locking, weight activated brake system.



COMES STANDARD

- Weight Activated Drag Brake (non-locking)
- Plastic Seat with built-in handle
- Front Wire Basket
- Flip-up Back
- Foldable
- 5" Wheels

OPTIONS

- Travel Bag
- Tray
- Rubber Seat Pad
- Basket Cover Bag
- Cup Holder
- Flashlight Attachment
- Cane Holder
- Designer Mobility Bags
- Colors: Blue, Red

WEIGHT CAPACITY: 250 lbs

WEIGHT: 12 lbs

SEAT HEIGHT: 21"

SEAT DIMENSIONS: 13.5"w x 13"d

HANDLE HEIGHT ADJUSTMENT: 31.5" - 36"

WIDTH BETWEEN HANDLES: 17.5"

OVERALL DIMENSIONS: 24.5"w x 22.5"d

APPROX. USER HEIGHT: 5'3" - 6'0"

BASKET DIMENSIONS: 9.5"l x 15"w x 6.5"d

DIMENSIONS WHEN FOLDED: 32"l x 24.5"w x 7"d



Cruiser Classic 4200C

The **Cruiser Classic** allows you to cruise along with a **weight activated** brake system. This walker **DOES NOT** have locking hand brakes for those unable to use them. Moderate-to-good stability is required to use the non-locking, weight activated brake system.



COMES STANDARD

- Weight Activated Drag Brake (non-locking)
- Padded Seat
- Under Seat Basket
- Flip-up Back
- Foldable
- 5" Wheels

OPTIONS

- Travel Bag
- Basket Cover Bag
- Cup Holder
- Flashlight Attachment
- Cane Holder
- Designer Mobility Bags
- Colors: Blue, Red

WEIGHT CAPACITY: 250 lbs

WEIGHT: 12 lbs

SEAT HEIGHT: 21"

SEAT DIMENSIONS: 11.5"w x 11.5"d

HANDLE HEIGHT ADJUSTMENT: 31.5" - 36.5"

WIDTH BETWEEN HANDLES: 17.5"

OVERALL DIMENSIONS: 24.5"w x 22.5"d

APPROX. USER HEIGHT: 5'3" - 6'0"

BASKET DIMENSIONS: 9"l x 15"w x 6.5"d

DIMENSIONS WHEN FOLDED: 32"l x 24.5"w x 8"d



Cruiser II 4201

The **Cruiser II** is ideal for beginning **gait training** with 5 inch front wheels and rear rubber tips to prevent the walker from rolling away. The rolling walker frame and seat provide the stability, safety, and confidence to **establish a natural walking pattern.**



Walker Skis
(Item # 40027GR)

COMES STANDARD

Large Plastic Seat
Front Wire Basket
Flip-up Back
Foldable
5" Wheels

OPTIONS

Travel Bag
Rubber Seat Pad
Cup Holder
Tray
Basket Cover Bag
Flashlight Attachment
Cane Holder
Designer Mobility Bags
Walker Skis
Color: Blue

WEIGHT CAPACITY: 250 lbs

WEIGHT: 10 lbs

SEAT HEIGHT: 20.25"

SEAT DIMENSIONS: 13.5"w x 13"d

HANDLE HEIGHT ADJUSTMENT: 32" - 35.5"

WIDTH BETWEEN HANDLES: 17.5"

OVERALL DIMENSIONS: 20.5"w x 20.5"d

APPROX. USER HEIGHT: 5'3" - 6'0"

BASKET DIMENSIONS: 9.5"l x 15"w x 6.5"d

DIMENSIONS WHEN FOLDED: 31.5"l x 20.5"w x 7"d



Cruiser De-Light 4010

The **Cruiser De-Light** is a real **pleasure!** It is the only 4-wheeled rolling walker that folds and has adjustable legs like a standard folding walker with the added benefits of locking hand brakes and a seat.



COMES STANDARD

Locking Hand Brakes
Plastic seat with built-in rubber seat pad
Foldable
5" Wheels

OPTIONS

Basket (Item #4010B)
Tray
Travel Bag
Seat Covers
Cup Holder
Flashlight Attachment
Cane Holder
Designer Mobility Bags
Colors: Blue, Red

WEIGHT CAPACITY: 250 lbs

WEIGHT: 16 lbs

SEAT HEIGHT: 19.5" - 23.5" **SEAT DIMENSIONS:** 15.5"w x 10.5"d

HANDLE HEIGHT ADJUSTMENT: 34" - 38"

WIDTH BETWEEN HANDLES: 18"

OVERALL DIMENSIONS: 27.5"w x 23"d

APPROX. USER HEIGHT: 5'4" - 6'2"

BASKET DIMENSIONS: 7"l x 20"w x 7.5"d

DIMENSIONS WHEN FOLDED: 35"l x 27.5"w x 9"d



Mini Mack 4214

Mack 4215

The **Mini Mack** is strong and reliable with a **400 pound weight capacity** and a lower seat height for more **petite** users.

The **Mack** is strong and reliable with a **400 pound weight capacity** and a standard height for **average height to tall** users.



COMES STANDARD

- Locking Hand Brakes
- Plastic seat with built-in rubber seat pad
- Front Wire Basket
- Flip-up Back
- Foldable
- 8" Wheels

OPTIONS

- Tray
- Basket Cover Bag
- Seat Cover
- Cup Holder
- Flashlight Attachment
- Cane Holder
- Designer Mobility Bags
- Colors: Blue, Red

COMES STANDARD

- Locking Hand Brakes
- Extra Wide Plastic Seat with built-in rubber seat pad
- Front Wire Basket
- Flip-up Back
- Foldable
- 8" Wheels

OPTIONS

- Tray
- Basket Cover Bag
- Seat Cover
- Cup Holder
- Flashlight Attachment
- Cane Holder
- Designer Mobility Bags
- Colors: Blue, Red

WEIGHT CAPACITY: 400 lbs **WEIGHT:** 20.5 lbs
SEAT HEIGHT: 19.5" **SEAT DIMENSIONS:** 16.75"w x 12.5"d
HANDLE HEIGHT ADJUSTMENT: 29.25" - 33.75"
WIDTH BETWEEN HANDLES: 20.5"
OVERALL DIMENSIONS: 24.5"w x 24.5"d
APPROX. USER HEIGHT: 4'11" - 5'5"
BASKET DIMENSIONS: 8"l x 17.5"w x 6"d
DIMENSIONS WHEN FOLDED: 35"l x 24.5"w x 11"d



WEIGHT CAPACITY: 400 lbs **WEIGHT:** 21.5 lbs
SEAT HEIGHT: 23" **SEAT DIMENSIONS:** 16.75"w x 12.5"d
HANDLE HEIGHT ADJUSTMENT: 34.5" - 39"
WIDTH BETWEEN HANDLES: 20.5"
OVERALL DIMENSIONS: 24.5"w x 27.5"d
APPROX. USER HEIGHT: 5'5" - 6'2"
BASKET DIMENSIONS: 10"l x 16.5"w x 6"d
DIMENSIONS WHEN FOLDED: 39"l x 24.5"w x 11"d



Mighty Mack 4216

The **strongest** of them all is the **Mighty Mack** with a **600 pound weight capacity** and a standard seat height for **average to tall** users.



COMES STANDARD

Locking Hand Brakes
Extra Wide Plastic Seat
with built-in rubber seat
pad
Front Wire Basket
Removable Back
Foldable
8" Wheels

OPTIONS

Tray
Basket Cover Bag
Cup Holder
Flashlight Attachment
Cane Holder
Designer Mobility Bags
Colors: Blue, Red

WEIGHT CAPACITY: 600 lbs

WEIGHT: 24 lbs

SEAT HEIGHT: 22"

SEAT DIMENSIONS: 20"w x 13.75"d

HANDLE HEIGHT ADJUSTMENT: 30.25" - 36.25"

WIDTH BETWEEN HANDLES: 23.5"

OVERALL DIMENSIONS: 27.5"w x 25.5"d

APPROX. USER HEIGHT: 5'5" - 6'2"

BASKET DIMENSIONS: 10"l x 16.5"w x 6"d

DIMENSIONS WHEN FOLDED: 28.5"l x 27.5"w x 11.5"d



NOVA Options

CLUTCHES

• Attach securely to walkers, wheelchairs, transport chairs, canes, and bed rails.

- Item # BGS01 (Black)
- Item # BGS02 (Blue)
- Item # BGS03 (Red)
- Item # BGS04 (Leopard)
- Item # BGS05 (Camouflage)



WRISTLETS

• Attach securely to walkers, wheelchairs, transport chairs, canes, and bed rails.

- Item # BGM01 (Black)
- Item # BGM03 (Pink)
- Item # BGM04 (Zebra)
- Item # BGM05 (Metropolitan)



Clutches and Wristlets are great for keys, cellphones, and other personal items.

HANGING WALKER POUCHES

- Fit all folding and 4200 series walkers
- Velcro to front of crossbar of folding walker or backrest of 4200 series
- Extra storage for lightweight objects



• Item # P4200SB (Aloha Blue)

• Item # 4001WP (Black)

MOBILITY SADDLE BAGS



• Attach securely to all walkers, wheelchairs, transport chairs, and bed rails

- Item # 4002AP (Aloha Pink)
- Item # 4002BF (Butterflies)
- Item # 4002BK (Black)
- Item # 4002CT (Cats)
- Item # 4002DG (Dogs)
- Item # 4002MP (Metropolitan)

Saddle Bags are so versatile! Use them with your mobility product or as an attractive over-the-shoulder bag.

NEW ACCESSORIES

NEW!

Glamour Line

ULTIMATE WALKER BAGS

• Versatile to fit on rolling walkers and folding walkers. Generous and roomy with convenient front pockets.

- Item # 4003BK (Black)
- Item # 4003HT (Houndstooth)
- Item # 4003EA (Edgy Aqua)
- Item # 4003LD (Leopard)



*be a
mobility
fashionista!*



NEW!

Glamour Line

MOBILITY HAND BAGS

• So sleek and stylish. Perfect to go from walker to over the shoulder. Comes with adjustable shoulder strap.

- Item # 4005BK (Black)
- Item # 4005EA (Edgy Aqua)
- Item # 4005HT (Houndstooth)
- Item # 4005LD (Leopard)



*be a
mobility
fashionista!*



NEW ACCESSORIES

NEW!

Glamour Line

UNIVERSAL MOBILITY BAGS

• This bag is global! - Fits on any bar, backrest, or handles. So many options, so many places to go. Large main compartment, multiple outer pockets.

- Item # 4004BK (Black)
- Item # 4004LD (Leopard)
- Item # 4004EA (Edgy Aqua)
- Item # 4004HT (Houndstooth)



*be a
mobility
fashionista!*



NEW!**Glamour Line****BACK AND SEAT COVERS**

- Add glam, add style, add luxurious feel. Make it all about you. Stand out from the crowd and make a statement!



• Item # 4007SL
(Snow Leopard)




• Item # 4007SC
(Safari Cheetah)

NEW!**CANE HOLDER**

- For rolling walkers and folding walkers
- Fits 1" diameter tubing and has multiple angle attachments
- Cane (with 7/8" diameter) snaps easily in and out of the holder
- Reversible mount for left or right hand application



• Item # CH-4000R



• Item # CH-1000



• Item # CH-2000



• Item # FL-1000



• Item # 4000TB



• Item # 4000T



- Easy-Access Velcro Top
- Water-Resistant



- Item # 4200RS (fits 4200/4201/4203/4208/4212)
- Item # 4230RS (fits 4200/4203/4212)

CUP HOLDERS

- Attach easily to walkers, wheelchairs, bed rails and transport chairs
- Adjusts to hold different sized items (CH-1000 only)
- Folds up when not in use (CH-1000 only)

FLASHLIGHT

- Adds extra light for safety in dark places
- Attaches to most walkers, transport chairs, and wheelchairs

TRAVEL BAG

- Fits all folding and rolling walkers
- Nylon bag, perfect for airline travel or storage
- 10"d x 36"w x 27"h

WALKER TRAY

- Carries food, drinks, and various personal belongings

BASKET COVER BAG

- Item # 4200P (for 4200/4201/4203/4212 basket)
- Item # 4208P (for 4208 basket)
- Item # 4218P (for 4202/4207 basket)
- Item # 4215P (for 4214/4215/4216 basket)

RUBBER SEAT PAD

- 1/2" thick rubber helps reduce back fatigue when resting
- Attaches to plastic seat with 4 screws

LIVE LIFE TO THE
FULLEST WITH YOUR NOVA:
The *Very Best* in Rolling Walkers.



Improve your overall mobility, safety, independence, and style by selecting from Nova's extensive line of rolling walkers and accessories.

WARNING



- Lock hand brakes before sitting and while seated
- Seat is for stationary sitting only. Do not ambulate or use as a transport chair while seated
- Do not use on stairs, escalators, or moving walkways
- Backrest must be attached to the walker prior to operating
- Consult your user manual before operating



as unique as you are

1-800-557-6682

www.NovaMedicalProducts.com

Authorized NOVA Dealer

