



NOVUM

Medical Products, Inc.

New Focus. New Commitment. New Innovative Solutions.

Bassinet Model NB-SSxB

Combines beauty with stability. Suited for all-around service in any nursery. Constructed of reinforced, heavy-gauge polished stainless steel for safety and durability.

Features:

- Stainless steel shelf is attached 8" from the ground
- NK Whisper-Glidesm 3" Casters
- Trendelenburg may be created from either end of the bassinet
- Extender shelf option available



Pictured with optional acrylic plastic basket and mattress pad



Dimensions: 31" L x 17.5" W x 37.75" H

Phone | 716.759.7200
Toll Free | 800.274.2742
Fax | 716.759.0700
Email | info@novummed.com
www.novummed.com



NOVUM

Medical Products, Inc.

New Focus. New Commitment. New Innovative Solutions.

Extender Shelf

Model 7028-EX



The stainless steel sliding extension shelf is available on any NK bassinet. It can be extended to a maximum of 31" to provide an infant changing surface. The shelf has 5/16" safety side rails supported by 1" tubular steel legs for added support.

Features:

- Shelf locks into place when extended
- Pull-out shelf measures 30" tall
- NK Whisper-Glidesm 3" Casters
- ★ Not supported on xDC models



Polycarbonate Basket

Model 69576X1-P

Clear polycarbonate with rigid one-piece construction. It is available for all NK bassinets.

Features:

- 3"x 5" card holder at head end
- Vent holes on bottom
- Dimensions: 26.5"L x 13"W 10"H



Bassinet Mattress

Model 7028X21

NK mattresses conform to all federal regulations. The outside cover is stain and odor resistant, anti-bacterial and flame retardant. Foam material meets California #117 code for flame retardant material. All mattresses are latex-free. Additional sizes available.

Dimensions: 25.75"L x 12"W x 2"H

Phone | 716.759.7200
Toll Free | 800.274.2742
Fax | 716.759.0700
Email | info@novummed.com
www.novummed.com