



System One
Sleep Therapy Systems

Category: CPAP

REMstar Pro



DS460S REMstar Pro C-Flex+
DS460HS REMstar Pro C-Flex+
with humidifier
DS460TS REMstar Pro C-Flex+
with Heated Tube
humidifier

Your home sleep testing and long-term CPAP solution

With **System One REMstar Pro** you now have the ability to assess your home sleep testing patients' CPAP prescription in Auto-Trial mode. **Auto-Trial mode** and **CPAP-Check mode** were designed to optimize pressure settings, while providing the long-term benefits of fixed CPAP therapy. At the same time, REMstar Pro includes virtually all of the features and programs that make System One the most prescribed sleep therapy system by US sleep centers.

System One features

- Advanced event detection
- C-Flex+ for comfort
- System One Resistance Control
- Mask fit check
- Waveform reports
- Oximetry module
- On-board memory
- Encore patient management system
- SD card and wired, wireless, and broadband modem capable

Humidification

- System One humidification control
- Pre-heat humidifier option
- Heated Tube humidifier option

Programs

- Flex Promise program
- Bi-level Rescue program
- System "One-Touch" program
- Partners in Compliance Management

REMstar Pro advantages

- CPAP, CPAP-Check, Auto-Trial modes
- C-Flex+ for comfort
- Advanced event detection
- Modem and SD card capable
- Optional, integrated oximetry available
- Heated Tube option

PHILIPS
RESPIRONICS

REMstar Pro specifications

Modes	Comfort features		Encore data capture and reporting			Reimbursement
	Flex comfort	System One Resistance Control	Patterns of use	Advanced event detection*	Flow waveform	
CPAP	C-Flex	Standard	Standard	Standard	Standard	E0601
CPAP-Check	C-Flex+					
Auto-Trial	A-Flex					

*Philips Respironics devices with advanced event detection detect, respond to (some models), and log these advanced events: Apnea Hypopnea Index (AHI), respiratory-related arousal (RERA), snore, flow limitation (FL), leak, vibratory snore (VS), obstructive airway apnea (OA), large leak (LL), clear airway apnea (CA), periodic breathing (PB), and hypopnea (H).

Specifications	
Pressure range	4-20 cm H ₂ O
Flex pressure relief	0 to 3
Ramp time	0 to 45 min (5-min increments)
Starting ramp pressure	4 to CPAP; patient adjustable (CPAP mode) 4 to min pressure (Auto mode)
Dimensions	7" L x 5.5" W x 4" H 18 cm L x 14 cm W x 10 cm H
Weight	3 lbs; 1.36 kg (without humidifier)
Filters	Foam and optional ultra-fine
Device set-up	LCD/control wheel/push button
Data storage capacity (minimum)	Display: compliance VIC, 1-, 7- and 30-day averages SD card > 1 year On-board: up to 6 months
Compliance meter	Breathing detection
Altitude compensation	Automatic
Electrical requirements	100-240 VAC, 50/60 Hz
Humidification	Integrates with System One 60 Series heated humidifier and System One 60 Series Heated Tube humidifier
Warranty	Two years (US only)
Miscellaneous	Lighted LEDs

Sleep therapy accessories	
Humidification	Part number
System One 60 Series humidifier	DS6H
System One 60 Series Heated Tube humidifier	DS6T
Water chamber (includes both 50 and 60 Series lids)	1063785
Heated Tube humidifier upgrade kit	6T-CONVKIT
Tubing	Part number
System One Heated Tube	SYS1HT15
System One Performance Tube (15mm)	1074400
Performance Tubing (22mm)	1032907
Filters	Part number
Pollen (foam) filter, reusable (1/pack)	1035443
Pollen (foam) filter, reusable (2/pack)	1029330
Ultra-fine filter, disposable (1/pack)	1063096
Ultra-fine filter, disposable (2/pack)	1035442
Ultra-fine filter, disposable (6/pack)	1029331
Accessories	Part number
Device side panel	1099592
Carrying case	1063857



Complies with the Buy American Act

Respironics, A-Flex, Auto-Trial, C-Flex, CPAP-Check, Encore, REMstar, and System One are trademarks of Koninklijke Philips Electronics N.V. and its affiliates. All rights reserved.

Please visit www.philips.us/remstarpro



© 2012 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

edoc JJ 10/30/12 MCI 4105269

Philips Healthcare is part of Royal Philips Electronics

Philips Respironics International Headquarters
+33 1 47 28 30 82
Philips Respironics Asia Pacific
+65 6882 5282
Philips Respironics Australia
+61 (2) 9947 0440
Philips Respironics United Kingdom
+44 800 1300 845

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)
www.philips.com/respironics