

Polyester Tipped Applicator

Plastic handle



- Small tip

- Inherent score mark

Actual Size

Description

This synthetic tipped applicator is made of medical-grade polyester fiber that is carded in our factory to maintain consistency. It is non-absorbent and contains no fatty acids or latent sugars, making it ideal for specimen transport, lab work, and diagnostic kits. It is bonded securely to a tapered plastic handle by an aqueous based adhesive. Inherent score mark on handle. Available sterile or non-sterile.

Dimensions

5.910" Overall Length (150.114 mm)

Tip

■ 0.078" Dia. (1.981 mm)

■ 0.687" Length (17.45 mm)

Handle

■ 0.085" Dia. tapered to .040" (2.159 to 1.016 mm)

■ 5.848" Length (148.539 mm)

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of compliance available with each shipment upon request

Qualities

- Medical grade quality
- Non-absorbent tip
- Contains no fatty acids or latent sugars
- Small tip for patient comfort
- Molded handle is semi-flexible

Applications

- Diagnostic specimen collection
- Recovery of micro-organisms and fluid
- Cell collection for DNA testing

Ref #	Packaging	
1000-PD <i>Non-sterile</i>	10/100	(10 bags per case, 100 swabs per bag)
25-1000 1PD <i>Sterile</i>	10/100/1	(10 boxes per case, 100 individually wrapped swabs per box, 1 per package)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com