

Foam Tipped Applicator

Polypropylene handle



- Soft, foam paddle tip

- Black polypropylene handle

Actual Size

Description

This soft, non-abrasive foam paddle tipped applicator is produced with reticulated 100ppi foam. It is heat staked to a black polypropylene handle.

Dimensions

6.044" Overall Length (153.518mm)

Tip

- 0.730" Length (18.542mm)
- 0.250" Width (6.35mm)

Handle

- 5.922" Length (150.419mm)
- 0.120" Dia. (3.048mm)

Qualities

- Lint free
- Good water absorption
- Soft non abrasive paddle tip
- No contaminating residues
- Can be used with IPA and DI water or dry

Standards

- ISO 9001:2008
- FDA GMP's

Applications

Industrial

- Applying adhesives or solutions
- Removing residues and particulates
- General cleaning

Ref #	Packaging
1276 BK	2/500 (2 poly bags per case, 500 per poly bag)



Specialty Applicators

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com