




# Puritan medical swabs – meeting healthcare needs for nearly a century.



## WHY PURITAN MEDICAL SWABS?

Puritan is a name known all over the world. We have built our **long-standing reputation** for **innovation** and **quality** by listening to our customers and meeting the evolving needs of the medical fields.











As always, Puritan responds to market needs with high-quality products. We have evolved alongside the healthcare field. Early in our history, the Puritan line was limited, primarily to tongue depressors and cotton tipped applicators. As medicine advanced, we kept pace, expanding our lines and even developing brand-new product categories and swabs for specialized medical applications.

The all-time Puritan best sellers shown above are time-tested examples. Look inside to see our standard product variety. And please remember, your feedback and custom orders are always welcome. The demands of the medical field are and will always be a driving force behind our R&D and our 100% commitment to consistent quality.



## COTTON TIPPED APPLICATORS

We are very particular about our raw materials. For these applicators, we only use US-grown cotton. We coil and card our USP-grade cotton fiber in our factory and inspect every swab to ensure optimum quality and consistency.

	Item No.	Item Description	Packaging pack, weight	Handle length x diameter	Tip length x diameter	Full Length
	803-PC	Polystyrene handle, standard tip.	10/10/100, 9 lbs	74.599 x 2.54 mm 2.937 x 0.1"	17.449 x 5.156 mm 0.687 x 0.203"	76.2 mm 3"
	25-803 2PC		10/100/2, 5 lbs			
	803-PC DBL	Polystyrene handle, double tipped.	10/10/50, 7 lbs	73.025 x 2.54 mm 2.875 x 0.1"	17.45 x 4.369 mm 0.687 x 0.172"	79.375 mm 3.125"
	803-PCL	Polystyrene handle, large tip.	10/5/50, 5 lbs	74.599 x 2.54 mm 2.937 x 0.1"	17.449 x 6.35 mm 0.687 x 0.25"	76.2 mm 3"
	803-WC	Wood handle, standard tip.	10/10/100, 10 lbs	73.812 x 2.54 mm 2.906 x 0.1"	15.875 x 4.775 mm 0.625 x 0.188"	76.2 mm 3"
	25-803 2WC		10/100/2, 5 lbs			
	803-WCL	Wood handle, large tip.	10/5/100, 7 lbs	73.812 x 2.54 mm 2.906 x 0.1"	17.449 x 5.943 mm 0.687 x 0.234"	76.2 mm 3"
	806-PC	Polystyrene handle, standard tip.	10/10/100, 17 lbs	147.624 x 2.54 mm 5.812 x 0.1"	17.449 x 5.156 mm 0.687 x 0.203"	152.4 mm 6"
	25-806 1PC		10/100/1, 6 lbs			
	25-806 2PC		10/100/2, 7 lbs			
	806-PCL	Polystyrene handle, large tip.	10/5/100, 10 lbs	149.225 x 2.54 mm 5.875 x 0.1"	17.449 x 6.35 mm 0.687 x 0.25"	152.4 mm 6"
	806-WC	Wood handle, standard tip.	10/10/100, 16 lbs	147.624 x 2.54 mm 5.812 x 0.1"	15.875 x 4.775 mm 0.625 x 0.188"	152.4 mm 6"
	25-806 1WC		10/100/1, 6 lbs			
	25-806 2WC		10/100/2, 7 lbs			
	25-806 5WC		10/20/5, 5 lbs			
	25-806 10WC		10/20/10, 6 lbs			
	806-WCL	Wood handle, large tip.	10/5/100, 10 lbs	147.624 x 2.54 mm 5.812 x 0.1"	17.449 x 5.943 mm 0.687 x 0.234"	152.4 mm 6"
	821-WC DBL	Wood handle, double tipped, tapered one end.	10/5/100, 9 lbs	147.624 x 2.108 mm 5.812 x 0.083"	15.875 x 4.775 mm 0.625 x 0.188"	152.4 mm 6"
	823-WC	Wood handle, tapered mini tip.	10/10/100, 8 lbs	73.812 x 2.108 mm 2.906 x 0.083"	12.7 x 3.048 mm 0.5 x 0.12"	76.2 mm 3"
	826-WC	Wood handle, tapered mini tip.	10/10/100, 12 lbs	147.624 x 2.108 mm 5.812 x 0.083"	12.7 x 3.048 mm 0.5 x 0.12"	152.4 mm 6"
	826-WC DBL	Wood handle, double tipped, tapered mini tip.	10/10/100, 13 lbs			

continued

## PACKAGING QUANTITY CODES




Our packaging codes work from the outside in. The first number represents the outermost package – this may be a case, a box, or other container. Next is what you'll find when you open the outer packaging – bags, wrappers, tubes, etc. Finally we get to the heart of the packaging – how many swabs are grouped together. Some examples follow:

10/50 = 10 bags of 50 swabs per bag, per case.

10/10/100 = 10 boxes of 10 bags per case, 100 swabs per bag.





10/100/1 = 10 boxes of 100 swabs with each swab individually wrapped or in a tube.

## COTTON TIPPED APPLICATORS, CONTINUED

	Item No.	Item Description	Packaging pack, weight	Handle length x diameter	Tip length x diameter	Full Length
	non-sterile   <b>sterile</b>					
	830-PC DBL	Paper handle, double tipped.	10/10/50, 8 lbs	73.025 x 2.54 mm 2.875 x 0.1"	17.449 x 4.368 mm 0.687 x 0.172"	76.2 mm 3"
	833-WCS	Wood handle, small tip, in plastic vials.	50/100, 8 lbs	73.812 x 2.108 mm 2.906 x 0.083"	15.875 x 3.962 mm 0.625 x 0.156"	76.2 mm 3"
	869-WCS	Wood handle, extra small tip.	10/10/100, 13 lbs	147.624 x 2.108 mm 5.812 x 0.083"	12.7 x 3.175 mm 0.5 x 0.125"	152.4 mm 6"

## OB/GYN AND PROCTOSCOPIC APPLICATORS








The best swabs – indeed, the best products in any category – are those that are designed to meet the real-world demands of the end user. These swabs are a good example.

	Item No.	Item Description	Packaging pack, weight	Handle length x diameter	Tip length x diameter	Full Length
	non-sterile   <b>sterile</b>					
	808	OB/GYN applicator, rayon tip, paper handle.	10/50, 5 lbs	203.2 x 3.962 mm 8 x 0.156"	28.575 x 12.7 mm 1.125 x 0.5"	207.975 mm 8.188"
	808 100		10/100, 8 lbs			
	808 BULK		500 bulk, 4 lbs			
	808 COTTON	OB/GYN applicator, cotton tip, paper handle.	10/50, 5 lbs	203.2 x 3.962 mm 8 x 0.156"	28.575 x 12.7 mm 1.125 x 0.5"	207.975 mm 8.188"
	25-808 1PR	OB/GYN applicator, rayon tip, polystyrene handle.	10/50/1, 8 lbs	196.85 x 3.81 mm 7.75 x 0.15"	28.575 x 12.7 mm 1.125 x 0.5"	203.2 mm 8"
	25-808 2PR		10/50/2, 11 lbs			
	816	Proctoscopic applicator, rayon tip, polypropylene handle.	10/50, 8 lbs	403.225 x 3.81 mm 15.875 x 0.15"	28.575 x 12.7 mm 1.125 x 0.5"	406.4 mm 16"
	816 100		10/100, 13 lbs			
	816 BULK		500 bulk, 6 lbs			

## CYTOLOGY COLLECTION DEVICES

Puritan cytology products are made of the very best materials, in the most effective shapes, to maximize sample collection and provide a comfortable patient experience at the same time.





Products below shown at 80% of actual size.

	Item No.	Item Description	Packaging pack, weight	Handle length x diameter	Tip length x diameter	Full Length
	2188	Histobrush, <sup>®</sup> soft nylon bristles, tapered tip.	10/100, 5 lbs	183.362 x 2.971 mm 7.219 x 0.117"	19.05 x 5.08 mm 0.75 x 0.2"	203.2 mm 8"
	25-2188		10/100/1, 8 lbs			
	2194	Rover's <sup>®</sup> spatula.	15/50, 4 lbs	167.5 x 4.4 mm 6.594 x 0.173"	30 x 16.7 mm 1.181 x 0.657"	196.85 mm 7.75"
	2195	Rover's <sup>®</sup> Cervex-Brush. <sup>®</sup>	35/25, 6 lbs	166.674 x 4.394 mm 6.562 x 0.173"	28.575 x 19.05 mm 1.25 x 0.75"	203.2 mm 8"
	2195 COMBI	Rover's <sup>®</sup> Cervex-Brush <sup>®</sup> combi.	4/25, 5 lbs	167.335 x 4.394 mm 6.588 x 0.173"	34.29 x 19.05 mm 1.35 x 0.75"	203.2 mm 8"
	2196	Rover's <sup>®</sup> EndoCervex-Brush. <sup>®</sup>	10/50, 3 lbs	167.487 x 4.394 mm 6.594 x 0.173"	33.782 x 5.435 mm 1.333 x 0.214"	203.2 mm 8"
	2197	CerviSoft, <sup>®</sup> soft foam tip.	10/100, 3 lbs	177.8 x 2.54 mm 7 x 0.1"	25.4 x 4.775 mm 1 x 0.188"	179.374 mm 7.062"
	25-2197		10/100/1, 7 lbs			
	2199	Histobrush, <sup>®</sup> soft nylon bristles, tapered tip.	10/100, 4 lbs	157.962 x 3.073 mm 6.219 x 0.121"	19.05 x 5.08 mm 0.75 x 0.2"	177.8 mm 7"
	25-2199		10/100/1, 8 lbs			

## CERVICAL SCRAPERS


Healthcare practitioners have counted on us for years to provide standard wooden scrapers. The quality of the wood, our manufacturing, and quality processes ensure a smooth finish for effective collection and patient comfort.

Products below shown at 30% of actual size.

	Item No.	Item Description	Packaging pack, weight	Full Length
	2173	Ayres style.	10/100, 3 lbs	177.8 mm 7"
	25-2173		10/100/1, 7 lbs	
	2181	Medical packaging style.	10/500, 19 lbs	177.8 mm 7"
	2183	Double-end.	10/500, 21 lbs	177.8 mm 7"
	2184	Lerner style.	10/500, 24 lbs	177.8 mm 7"





## FOAM WOUND MEASUREMENT DEVICES

The soft foam-tipped applicator is flexible to conform to wound shape and features a metric rule for easy measuring. Perfect for forensics, emergency rooms, long-term care, and more. Choose sterile or non-sterile.

	Item No.	Item Description	Packaging pack, weight	Handle length x diameter	Tip length x diameter	Full Length
	1506-PFB DM	DM stick, polystyrene handle, foam tip.	20/50, 2 lbs	149.225 x 2.54 mm 5.875 x 0.1"	9.98 x 5.08 mm 0.393 x 0.2"	152.4 mm 6"
	25-1506 1PF DM		4/50/1, 2 lbs			

## APPLICATOR STICKS

Made for general purpose functionality. Choose the size and shape that works best for your applications – specimen transfer, mixing rods, and laboratory work. Made with straight cut ends.

	Item No.	Item Description	Packaging pack, weight	Dimensions length x diameter
	805	Wood.	12/1000, 23 lbs	304.8 x 2.108 mm 12 x 0.083"
	25-807 2W		20/1000, 18 lbs	147.624 x 2.108 mm 5.812 x 0.083"
	807	Wood.	10/100/2, 6 lbs	304.8 x 3.175 mm 12 x 0.125"
	809		12/500, 25 lbs	
	811	Wood.	2/1000, 12 lbs	457.2 x 3.175 mm 18 x 0.125"

## CUTICLE/ORANGE STICKS

Single and double bevel ends offer numerous choices for all kinds of applications. Excellent for podiatric practices and long-term care. Ideal as an all-purpose stick for most clinical, diagnostic, and forensic settings. Products below shown at 40% of actual size.

	Item No.	Item Description	Packaging pack, weight	Dimensions length x diameter
	2900	Wood, double bevel both ends.	5000 bulk, 16 lbs	177.8 x 3.962 mm 7 x 0.156"
	2910		10/100, 5 lbs	
	2911	Wood, double bevel both ends, wrapped.	4/125/1, 6 lbs	177.8 x 3.962 mm 7 x 0.156"
	2920	Wood, single bevel, double bevel ends.	5000 bulk, 18 lbs	
	2925	Wood, single bevel, pointed one end.	10/100, 4 lbs	152.4 x 3.962 mm 6 x 0.156"
	2930	Wood, single bevel, double bevel ends.	10/100, 5 lbs	177.8 x 3.962 mm 7 x 0.156"
	2940	Wood, double bevel both ends.	2500 bulk, 20 lbs	177.8 x 4.775 mm 7 x 0.188"
	2950		10/100, 6 lbs	
	2970	Wood, single bevel, double bevel ends.	10/100, 6 lbs	177.8 x 4.775 mm 7 x 0.188"

## TONGUE DEPRESSORS

While these are still called tongue depressors, wooden blades are used for a wide variety of applications, from patient oral exams to mixing or applying medicines, to crafts and industrial use.



Products below shown at 40% of actual size.

	Item No.	Item Description	Packaging pack, weight	Dimensions length x width
	non-sterile   sterile			
	700	Infant.	10/1000, 37 lbs	114.3 x 9.53 mm 4.5 x 0.375"
	25-700		10/100/1, 8 lbs	
	704	Standard.	10/500, 35 lbs	152.4 x 17.475 mm 6 x 0.688"
	704 BULK		3,700 bulk, 28 lbs	
	705	Regular.	10/500, 38 lbs	152.4 x 19.05 mm 6 x 0.75"
	25-705		10/100/1, 10 lbs	
	706	Standard.	12/100, 9 lbs	152.4 x 17.475 mm 6 x 0.688"
	709	Standard, individually wrapped.	10/250, 21 lbs	152.4 x 17.475 mm 6 x 0.688"
	710	Junior, individually wrapped, plain or flavored.	10/250, 21 lbs	139.7 x 15.875 mm 5.5 x 0.625"
	710 CHERRY		20/50, 10 lbs	
	710 GRAPE			
	710 STRAWBERRY			
	711	Junior.	10/500, 29 lbs	139.7 x 15.875 mm 5.5 x 0.625"
	25-711		10/100/1, 10 lbs	
	712		12/100, 8 lbs	

## MEDICAL SPOONS

Today Puritan Medical Products manufactures some of the most sophisticated single-use medical and diagnostic products on the market. But we still pay close attention to our simple wood products.

Products below shown at 70% of actual size.

	Item No.	Item Description	Packaging pack, weight	Dimensions length x diameter
	625	Wood spoon, double end.	10/1000, 32 lbs	76.2 x 19.456 mm 3 x 0.766"
	645	Wood spoon, double end.	10/500, 31 lbs	127 x 25.4 mm 5 x 1"

## ECONOMY GRADE HOSPITAL & BRIGHTWOOD LINES

These lines offer lower-cost alternatives and include both US-made and imported materials. Spoons, blades, swabs, and applicators are ideal for mixing and dispensing medicines – handy for clinical kits and field work. Products below shown at 50% of actual size.

	Item No.	Item Description	Packaging pack, weight	Dimensions length x width
	non-sterile   <b>sterile</b>			
	625 HOSPITAL	Wood spoon, double end. Hospital® brand.	10/1000, 34 lbs	76.2 x 19.456 mm 3 x 0.766"
	704 BRIGHTWOOD	Standard size tongue blade. Brightwood® brand.	10/500, 35 lbs	152.4 x 17.475 mm 6 x 0.688"
	704 HOSPITAL	Standard size tongue blade. Hospital® brand.	10/500, 35 lbs	152.4 x 17.475 mm 6 x 0.688"
	25-705 HOSPITAL		10/100/1, 11 lbs	
	711 HOSPITAL	Junior size tongue blade. Hospital® brand.	10/500, 29 lbs	139.7 x 15.875 mm 5.5 x 0.625"
	25-711 HOSPITAL		10/100/1, 9 lbs	

	Item No.	Item Description	Packaging pack, weight	Handle length x diameter	Tip length x diameter	Full Length
	non-sterile   <b>sterile</b>					
	803-WC HOSPITAL	Standard cotton tip, wood handle. Hospital® brand.	10/10/100, 8 lbs	73.812 x 2.108 mm 2.906 x 0.083"	15.875 x 4.775 mm 0.625 x 0.188"	76.2 mm 3"
	25-803 2WC HOSPITAL		10/100/1, 11 lbs			
	806-WC HOSPITAL	Standard cotton tip, wood handle. Hospital® brand.	10/10/100, 8 lbs	147.624 x 2.108 mm 5.812 x 0.083"	15.875 x 4.775 mm 0.625 x 0.188"	152.4 mm 6"
	25-806 1WC HOSPITAL		10/100/1, 6 lbs			
	25-806 2WC HOSPITAL		10/100/2, 7 lbs			




  

	Item No.	Item Description	Packaging pack, weight	Dimensions length x diameter
	807 HOSPITAL	Wood applicator. Hospital® brand.	20/1000, 19 lbs	147.624 x 2.108 mm 5.812 x 0.083"

## POPULES®

Super-convenient pre-filled popules with self-saturating, medical grade foam heads. Specify one of our standard solutions or let us know what you need.

Products below shown at 60% of actual size.

	Item No.	Item Description	Packaging pack, weight	Handle length x diameter	Tip length x diameter	Full Length
	4100	Self-saturating foam swab, 3% hydrogen peroxide.	4/3/50, 6 lbs	84.124 x 6.858 mm 3.312 x 0.27"	28.575 x 15.875 mm 1.250 x 0.625"	90.474 mm 3.562"
	4600 P	Self-saturating foam swab, 10% povidone iodine prep solution.	4/150/1, 7 lbs	84.124 x 6.858 mm 3.312 x 0.27"	34.925 x 17.449 mm 1.375 x 0.687"	90.474 mm 3.562"
	4620	Self-saturating foam swab, 70% isopropyl alcohol and deionized water.	4/150/1, 7 lbs	84.124 x 6.858 mm 3.312 x 0.27"	34.925 x 17.449 mm 1.375 x 0.687"	90.474 mm 3.562"





**Puritan Products**  
A Green Company

**PURITAN MEDICAL PRODUCTS CO. LLC**  
P.O. BOX 149, 31 SCHOOL STREET  
GUILFORD, MAINE 04443-0149

Tel 800-321-2313 (US and Canada)  
207-876-3311

Fax 800-323-4153 (US and Canada)  
207-876-3130

Web [www.puritanmedproducts.com](http://www.puritanmedproducts.com)  
Email [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

#### **PURITAN IS PASSIONATE ABOUT BEING GREEN**

We think green and give energy efficiency and resource stewardship high priority. Below are just a few things that make us green. Visit our web site for more.

- Puritan runs on 100% renewable energy (wood process waste) for all building heat – all building hot water – and all manufacturing-process steam conditioning and drying.
- We have an ongoing energy upgrade initiative, replacing fixtures and equipment with more efficient solutions.
- Puritan packaging materials are made up of at least 65% post-consumer recycled content.

#### **THE PURITAN ADVANTAGE™**

Every one of our products in the Puritan and PurSwab lines comes with the Puritan Advantage™ built right in.

We are a US manufacturer, with R&D, fabrication, sterilization, packaging, shipping, QC and oversight, sales, and customer service under one roof. Our products are consistently reliable and readily available.

We are a global company. Our medical devices carry the CE mark.

ISO 9001:2008 ISO 13485:2003 CE