

Foam Tipped Applicator

Polypropylene shaft



- Rectangular foam tip

- Rigid rectangular polypropylene shaft

Description

This foam tipped applicator is produced with reticulated 100ppi Medical Grade foam. It is thermally bonded to a rigid rectangular polypropylene shaft.

Dimensions

5.188" Overall Length (131.78mm)

Tip

- 1" Length (25.4mm)
- .531" Width (13.49mm)

Shaft

- 5" Length (127mm)
- .125" Dia. (3.18mm)
- .195" Width (4.95mm)

Qualities

- Lint and residue free
- Soft and non-abrasive
- Can be used with IPA and DI water, or dry
- Large tip

Standards

- ISO 9001:2008
- FDA GMP's

Applications

Industrial

- Buccal collection
- Applying solutions or solvents
- General cleaning

Ref #	Packaging
1605-PSF RECT	20/50 (20 resealable poly bags per case, 50 swabs per bag)
1605-PSFB RECT	1,000 bulk (1 poly bag per case, 1000 swabs per bag)



Puritan®
Quality since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com