



Puritan
Quality since 1919

800-321-2313
PURITANMEDPRODUCTS.COM

DATA SHEET

Foam Tipped Applicator

Large round foam tip, polypropylene handle.

This foam tipped applicator is made with 100PPI medical grade foam. The round head offers a large surface area for effective collection of small samples. This applicator is ideal for a variety of clinical and DNA applications.

Dimensions

5.125" Overall Length (130.175 mm)

Handle

- 0.125" Dia. (3.175 mm)
- 5.0" Length (127.0 mm)

Tip

- 0.906" Dia. (23.012 mm)
- 1.125" Length (28.575 mm)

Qualities

- Reticulated 100PPI foam tip
- Round lollipop shape
- Good water absorption
- Lint and residue free
- Puritan's PurSwab brand
- Made in USA

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Applications

- Removing residues or particulates
- Can be used with IPA and DI water, or dry
- General cleaning swab



Ref#	Packaging
1805-PF RND (Non-sterile)	20/50 (20 bags, 50 swabs per bag)

