



**Puritan**  
Quality since 1919

800-321-2313  
PURITANMEDPRODUCTS.COM

#### DATA SHEET

## Foam Tipped Applicator

Foam-over-cotton tip, polypropylene handle.

This absorbent foam tipped applicator is produced with reticulated 100PPI medical grade foam over high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is thermally bonded to a polypropylene handle ensuring the sponge-like foam tip cleans up unwanted particles without releasing any into your critical environment.

### Dimensions

6.125" Overall Length (155.575 mm)

#### Handle

- 0.1" Dia. (2.54 mm)
- 6.00" Length (152.40 mm)

#### Tip

- 0.375" Dia. (9.525 mm)
- 0.937" Length (23.800 mm)

### Qualities

- Lint and residue free
- Cotton under foam provides extra absorbency
- Can be used with IPA and DI water, or dry
- Made in USA

### Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

### Applications

- General cleaning swab
- Removing contamination
- Applying adhesives and other solutions



Ref#	Packaging
1806-PCF (Non-sterile)	20/50 (20 bags per box, 50 swabs per bag)



**Puritan**