



Puritan
Quality since 1919

800-321-2313

PURITANMEDPRODUCTS.COM

DATA SHEET

Foam Tipped Applicator

Standard foam tip, polypropylene handle.

This foam tipped applicator is made with reticulated 100PPI medical grade foam. It is thermally bonded to a polypropylene handle ensuring the sponge-like foam tip cleans up unwanted particles without releasing any into your critical environment.

Dimensions

6" Overall Length (152.4 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Tip

- 0.188" Dia. (4.775 mm)
- 0.625" Length (15.875 mm)

Qualities

- Reticulated 100PPI medical grade foam tip
- Soft, non-abrasive tip
- Can be used dry or with IPA and DI water
- Lint and residue free
- Made in USA

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Applications

- Critical environment industry
- Removing residues and particles
- Applying adhesives or other solutions
- General cleaning



Ref#	Packaging
1806-PF (Non-sterile)	20/50 (20 bags per box, 50 swabs per bag)



Puritan®