

1806-WCFL PurSwab Foam Tipped Applicator



Features and Benefits

Tip: <ul style="list-style-type: none"> • Reticulated 100PPI medical grade foam • 25.38mm length • 9.52mm width • Large foam over cotton tip 	Shaft: <ul style="list-style-type: none"> • Northern white birch • 152.28mm length • 2.54mm diameter 	Bond: <ul style="list-style-type: none"> • Friction 	Uses: <ul style="list-style-type: none"> • General cleaning swab 	Packaging: <ul style="list-style-type: none"> • 1,000 per case 20 bags of 50 in recloseable poly zip bags • Also available in bulk pack 4 poly bags of 250 Item: 1806-WCFLB 	Qualities: <ul style="list-style-type: none"> • Lint and residue free • Cotton under foam provides extra absorbancy • Soft
---	--	---	--	--	--



Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com