

# Foam Tipped Applicator

Wood handle



- Chemical resistant foam over cotton tip

- White Birch Shaft

Actual Size

## Description

This slightly abrasive foam tipped applicator is produced with chemical resistant foam over high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a rigid white birch handle.

## Dimensions

**6.062" Overall Length (153.975 mm)**

### Tip

- 0.812" Length (20.625 mm)
- 0.375" Width (9.525 mm)
- 0.3" Thickness (7.62 mm)

### Handle

- 5.812" Length (147.625 mm)
- 0.100" Dia. (2.54 mm)

## Standards

- ISO 9001:2008
- FDA GMP's

## Qualities

- Slightly textured foam for light abrasion
- Cotton under foam provides extra absorbency
- Will not swell or break apart with most solvents
- Shaft made from a renewable resource

## Applications

### Critical Environment / Industrial

- Applying adhesives or other solutions
- Removing dirt, contamination, fluxes or resins
- General cleaning

Ref #	Packaging
<b>1806-WCF CR</b>	20/50 (20 resealable poly bags per case, 50 swabs per bag)



Specialty Applicators

Puritan Medical Products Company LLC  
Guilford, Maine 04443-0149 USA

**Tel** 800-321-2313 (US and Canada Only)  
207-876-3311

**Fax** 800-323-4153  
207-876-3130

**Web** [www.puritanmedproducts.com](http://www.puritanmedproducts.com)

**Email** [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

**1-800-321-2313**  
[puritanmedproducts.com](http://puritanmedproducts.com)