



Puritan
Quality since 1919

800-321-2313

PURITANMEDPRODUCTS.COM

DATA SHEET

Foam Tipped Applicator - Dry Transport System

Acetal handle - sterile

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kits. This speciality tipped applicator is made with 100ppi medical grade polyurethane foam. It is boned securely to a tapered acetal handle.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.062" (152.4 mm)

Swab Handle

- 0.098" Dia. (2.489 mm) tapered to 0.049" (1.321 mm)
- 5.969" Length (152.298 mm)

Swab Tip

- 0.141" Dia. (3.581 mm)
- 0.547" Length (13.894 mm)

Qualities

- Foam has excellent collection and release properties
- Soft, small tip for patient comfort
- Lint and residue free

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Applications

- Specimen collection
- Diagnostic sampling
- Recovery of microorganisms and fluids



Ref#	Packaging	
25-1406 1PF BT	5/100/1	5 boxes per case, 100 tubes per box, 1 swab per tube