


25-1607 1PF SC PurCollect Environmental Sampling Collection Swab



Features and Benefits

Tip: <ul style="list-style-type: none"> • Open cell Macrofoam • Large surface area • 0.438" Width • 1" Length • Thermally bonded 	Shaft: <ul style="list-style-type: none"> • White polypropylene scored at 1.59" from tip for easy head detachment • Prominent thumb to indicate hand placement • Ensures uniform pressure applied 	Standards: <ul style="list-style-type: none"> • ISO 13485:2003 • FDA GMP's • Certificate of Sterility and/or Compliance available with each shipment upon request 	Uses: <ul style="list-style-type: none"> • Environmental sample collection for microbial detection and decontamination 	Packaging: <ul style="list-style-type: none"> • 10 boxes of 50 individually wrapped packages per case, 1 per package (10/50/1) • Sterile 	Qualities: <ul style="list-style-type: none"> • Sterile • Single pack to avoid cross contamination of swab stock • Foam has excellent collection and release properties • Validated in swab sampling protocol for the detection of bacillus anthracis
--	---	---	--	---	--



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com