



Puritan
Quality since 1919

800-321-2313

PURITANMEDPRODUCTS.COM

DATA SHEET

HydraFlock® Flocked Swab - Dry Transport System

Standard tip, polystyrene handle with breakpoint at 80mm - sterile

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kits. The patented HydraFlock® flocked swab provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.04" (153.416 mm)

Swab Handle

- 0.1" Dia. (2.54 mm) tapered to 0.066" (1.676 mm)
- 6.01" Length (152.654 mm)

Swab Tip

- 0.2" Dia. (5.08 mm)
- 0.69" Length (17.526 mm)

Qualities

- Medical grade quality
- Highly absorbent capacity
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA

Standards

- CE Marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Applications

- For collection and release of biological specimens
- Buccal cell collection
- Specimen collection
- Forensic evidence collection
- Rapid diagnostic collection



Ref#	Packaging	
25-3306-H BT	5/100/1	5 boxes per case, 100 tubes per box, 1 swab per tube