



Puritan
Quality since 1919

800-321-2313

PURITANMEDPRODUCTS.COM

DATA SHEET

PurFlock Ultra® Flocked Swab - Dry Transport System

Standard tip, polystyrene handle with breakpoint at 80mm - sterile

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kits. PurFlock Ultra® flocked swabs provide excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.04" (153.416 mm)

Swab Handle

- 0.1" Dia. (2.54 mm) tapered to 0.066" (1.676 mm)
- 6.01" Length (152.654 mm)

Swab Tip

- 0.2" Dia. (5.08 mm)
- 0.69" Length (17.526 mm)

Qualities

- Medical grade quality
- Highly absorbent tip
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA

Standards

- CE Marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Applications

- For collection of bacteria or microorganisms
- Buccal cell collection
- Specimen collection
- Forensic evidence collection
- Rapid diagnostic testing
- DNA Testing



Ref#	Packaging	
25-3306-U BT	5/100/1	5 boxes per case, 100 tubes per box, 1 swab per tube