

PurFlock® Ultra Flocked Swab

Polystyrene handle, mini-tip



- Mini flock tip

- Polystyrene handle

Actual Size



**PurFlock®
Ultra**

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Description

An ideal choice in flock swabs with excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests. Patents pending: Europe 10251389.2, US 12/849,250.

Dimensions

6.04" Overall Length (153.416 mm)

Tip

■ 0.135" Dia. (3.429 mm)

■ 0.725" Length (18.415 mm)

Handle

■ 0.098" tapered to 0.049" Dia.
(2.489 tapered to 1.245 mm)

■ 6" Length (152.4 mm)

Standards

■ CE marked for medical use

■ ISO 9001:2008

■ ISO 13485:2003

■ FDA GMP's

■ Certificate of sterility and/or conformance available with each shipment upon request

Qualities

■ Medical grade quality

■ Highly absorbent tip

■ Soft smooth tip

Applications

■ For Nasopharyngeal and Urethrogenital cell collection

■ Rapid diagnostic testing

■ Specimen collection

■ Forensic evidence collection

■ For collection of bacterial or micro-organisms

Ref #	Packaging
3316-U <i>Non-sterile</i>	10/100 (10 bags per case, 100 per bag)
25-3316-U <i>Sterile</i>	10/50/1 (10 boxes per case, 50 individually wrapped packages, 1 swab per package)
25-3316-U TT MC <i>Sterile</i>	10/50/1 (10 boxes per case, 50 individually wrapped transport devices, 1 swab per device)

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com