


PurFlock® Ultra Flocked Swab

Polystyrene handle with molded break point, micro ultrafine tip



- Micro ultrafine flock tip

- Polystyrene handle with molded break point

Actual Size



**PurFlock®
Ultra**

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Description

An ideal choice in flock swabs with excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests. Patents pending: Europe 10251389.2, US 12/849,250.

Dimensions

6.03" Overall Length (153.162 mm)

Tip

■ 0.140" Dia. (3.356 mm)

■ 0.312" Length (7.925 mm)

Handle

■ 0.100 tapered to 0.036" Dia.
(2.54 tapered to 0.914 mm)

■ 5.99" Length (152.146 mm)

Standards

■ CE marked for medical use

■ ISO 9001:2008

■ ISO 13485:2003

■ FDA GMP's

■ Certificate of sterility and/or conformance available with each shipment upon request

Qualities

■ Medical grade quality

■ Highly absorbent tip

■ Soft smooth tip

Applications

■ For Nasopharyngeal and Urethral cell collection

■ Rapid diagnostic testing

■ Specimen collection

■ Forensic evidence collection

Ref #	Packaging
3318-U <i>Non-sterile</i>	10/100 (10 bags per case, 100 per bag)
25-3318-U <i>Sterile</i>	10/50/1 (10 boxes per case, 50 individually wrapped packages, 1 swab per package)
25-3318-U TT MC <i>Sterile</i>	10/50/1 (10 boxes per case, 50 individually wrapped transport devices, 1 swab per device)

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com