

HydraFlock® Flocked Swabs

Polystyrene handle with molded break point, large tip



- Large flock tip

- Polystyrene handle with molded break point

Actual Size



HydraFlock®

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Description

A new era in flock swab technology. HydraFlock provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens. Patents pending: Europe 102513888.4, US 12/849,343.

Dimensions

6.06" Overall Length (153.924 mm)

Tip

■ 0.212" Dia. (5.385 mm)

■ 0.640" Length (16.256 mm)

Handle

■ 0.100 tapered to 0.066" Dia.
(2.54 tapered to 1.676 mm)

■ 6.01" Length (152.654 mm)

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Qualities

- Medical grade quality
- Highly absorbent tip
- Large, soft smooth tip

Applications

- For Buccal cell collection
- Rapid diagnostic testing
- Specimen collection
- Forensic evidence collection
- DNA testing

| Ref # | Packaging |
|--|--|
| 3406-H <i>Non-sterile</i> | 10/100 (10 bags per case, 100 per bag) |
| 25-3406-H <i>Sterile</i> | 10/50/1 (10 boxes per case, 50 individually wrapped packages, 1 swab per package) |
| 25-3406-H TT MC <i>Sterile</i> | 10/50/1 (10 boxes per case, 50 individually wrapped transport devices, 1 swab per device) |

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com