


# Polyester Tipped Applicator

Polystyrene handle



- Mini tip

- Ultrafine polystyrene handle

Actual size

## Description

This specialty applicator is tipped with medical grade polyester that we card ourselves. They're non-absorbent and contain no fatty acids or latent sugars, making them ideal for transport, lab work, and diagnostic kits. The Miniature tip is ideal for male specimen collection including STD testing.

## Dimensions

6.062" Overall Length (153.975mm)

### Tip

■ 0.078" Dia. (1.981mm)

■ 0.687" Length (17.45mm)

### Shaft

■ 0.098" / 0.052" Dia. (2.489mm / 1.321mm)

■ 4" / 2" Length (101.6mm / 50.8mm)

## Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or compliance available with each shipment upon request

## Qualities

- Medical grade quality
- Synthetic polyester contains no fatty acids or latent sugars
- Thin flexible polystyrene handle allows for minimal patient discomfort

## Applications

- Male specimen collection
- STD testing
- Testing and screening swab

Ref #	Packaging
<b>25-800 1PD 50</b> (sterile)	10/50/1 (10 boxes of 50 individually wrapped packages; 1 per package)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC  
Guilford, Maine 04443-0149 USA

**Tel** 800-321-2313 (US and Canada Only)  
207-876-3311

**Fax** 800-323-4153  
207-876-3130

**Web** [www.puritanmedproducts.com](http://www.puritanmedproducts.com)

**Email** [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

**1-800-321-2313**  
**[puritanmedproducts.com](http://puritanmedproducts.com)**